

IC Application Note No. C-113

Title:	Determination of lysine and standard		
	cations in a lysine sample		

Summary: Determination of lysine as well as sodium, ammonium, potassium and calcium in bulk lysine using cation chromatography with direct conductivity detection.


Sample:Bulk lysineSample Preparation:50 mg diluted in 100 mL ultrapure water

Column: 6.1010.220 Metrosep C 2 – 150

Eluent: 2.7 mmol/L oxalic acid

Flow: 1.0 mL/min

Injection Volume: 20 µL



Results:	Sodium	Ammonium	Potassium	Calcium	Lysine
	mg/kg	mg/kg	mg/kg	mg/kg	%
	2.1	0.14	0.07	0.04	98