
IC Application Note CIC–024

Organic chloride in naphtha fraction of crude oil distillation according to ASTM D8150



The content of organic chloride in crude oil is determined according to ASTM D8150 in the naphtha fraction after distillation. The naphtha fraction is whashed with caustic and water, respectively, to remove hydrogen sulfide and inorganic halides. Here, the determination of organic chloride after inline combustion is presented. Although the sulfur content was of no interest in this application, the same setup allows sulfur quantification.

Results

| | Mean [mg/kg] (N = 3) | RSD [%] (N = 3) |
|------------|-------------------------|--------------------|
| 1 Chlorine | 0.4 | 0.7 |
| 2 Sulfur | 146.3 | 0.5 |



Sample

Crude oil distillate

Sample preparation

The sample is analyzed by Combustion IC with flame sensor technology and intelligent Partial Loop Injection Technique with Inline Matrix Elimination.

Columns

| Metrosep A Supp 5 - 150/4.0 | 6.1006.520 |
|-----------------------------|------------|
| Metrosep A Supp 5 Guard/4.0 | 6.1006.500 |
| Metrosep A PCC 2 HC/4.0 | 6.1006.340 |

Solutions

| Eluent | 3.2 mmol/L sodium carbonate 1.0 mmol/L sodium hydrogen carbonate |
|-----------------------|---|
| Suppressor regenerant | 100 mmol/L sulfuric acid |
| Rinsing solution | STREAM |
| Absorber solution | 100 mg/L hydrogen peroxide |

Parameters

| Flow rate | 0.7 mL/min |
|-----------------------|--------------|
| Injection volume (IC) | 50 µL (MiPT) |
| P _{max} | 15 MPa |
| Recording time | 20 min |
| Column temperature | 30 °C |

Combustion parameters

| Argon | 100 mL/min |
|--|------------|
| Oxygen | 300 mL/min |
| Oven temperature | 1050 °C |
| Post-combustion time | 120 s |
| Initial volume of absorption solution | 2.0 mL |
| Water inlet | 0.1 mL/min |
| | |


Analysis


Conductivity after sequential suppression

Instrumentation

| 930 Compact IC Flex Oven/SeS/PP/Deg | 2.930.2560* |
|--|-------------|
| IC Conductivity Detector | 2.850.9010* |
| MSM Rotor A | 6.2832.000* |
| Adapter sleeve for Suppressor Vario | 6.2842.020* |
| 920 Absorber Module | 2.920.0010* |
| Combustion Module (oven and ABD) | 2.136.0700* |
| Autosampler MMS 5000 | 2.136.0800 |
| Kit for liquid sampling | 6.7303.000 |

* available as 930 Metrohm Combustion IC (2.930.9010)





www.metrohm.com

