
## IC Application Note CIC-028

# Fluorine and chlorine in iron ore by Combustion Ion Chromatography



Iron ore is an important resource for steel production. Its natural content of halogens is a quality characteristic due to the corrosiveness of the respective halogenides. Combustion IC applying the sacrificial vial technology is used for the analysis of fluorine and chlorine in ore.  $WO_3$  usually is added to improve the release of  $SO_2$  and therefore sulfur recovery. In this application, it also significantly improves the recovery of fluoride.

### Results

	Concentration [%] (N = 3)	RSD [%] (N = 3)	Recovery CRM [%] (no WO <sub>3</sub> )	Recovery CRM [%] (Sample : $WO_3 = 1:1$ )
1 Fluorine	84.0	1.7	80.0	92.7
2 Chlorine	70.2	1.1	94.7	97.8

Peaks 3 and 4 correspond to bromine and sulfur (not quantified)



#### Sample

Iron ore.

#### **Sample preparation**

The sample (5 to 15 mg) and in minimum the same amount of  $WO_3$  is placed in a quartz tube with quartz wool and subsequently fixed in the vial with quartz wool. This sacrificing vial is analyzed by Combustion IC with flame sensor technology and intelligent Partial Loop Injection Technique with Inline Matrix Elimination.

#### **Columns**

Metrosep A Supp 16 - 150/4.0	6.1031.420
Metrosep A Supp 16 Guard/4.0	6.1031.500
Metrosep A PCC 2 HC/4.0	6.1006.340
Metrosep I Trap 1 - 100/4.0	6.1014.200
Metrosep A Trap 1 - 100/4.0	6.1014.000

#### **Solutions CIC**

<u>Eluent</u>	7.5 mmol/L sodium carbonate 0.75 mmol/L sodium hydroxide	
Suppressor regenerant	100 mmol/L sulfuric acid	
Rinsing solution	STREAM	
Absorber solution	100 mg/L H <sub>2</sub> O <sub>2</sub>	

#### **Analysis**

Conductivity after sequential suppression

#### **Parameters**

Flow rate	0.8 mL/min
Injection volume (IC)	200 μL (MiPT)
P <sub>max</sub>	20 MPa
Recording time	18 min
Column temperature	45 °C


#### **Combustion parameters**


Argon	100 mL/min	
Oxygen	300 mL/min	
Oven temperature	1100 °C	
Post-combustion time	500 s	
Initial volume of absorption solution	2.0 mL	
Absorber solution feed	0.2 mL/min	
Water inlet	0.2 mL/min	
Post-combustion rinsing volume	1.0 mL	

#### Instrumentation

930 Compact IC Flex Oven/SeS/PP/Deg	2.930.2560*
IC Conductivity Detector	2.850.9010*
MSM Rotor A	6.2832.000*
Adapter sleeve for Suppressor Vario	6.2842.020*
920 Absorber Module	2.920.0010*
Combustion Module (oven and ABD)	2.136.0700*
Autosampler MMS 5000	2.136.0800
Kit for solid sampling	6.7302.000

<sup>\*</sup> available as 930 Metrohm Combustion IC (2.930.9010)





www.metrohm.com

