Thermo. Titr. Application Note No. H-003

Title:	Determination of Sulfate in Phosphoric Acid	
r		
Scope:	Determination of the sulfate content of wet process phosphoric acid.	
Principle:	An aliquot of concentrated phosphoric acid is titrated with standard barium chloride solution to a single thermometric endpoint.	
Reagent:	Standard c (BaCl ₂) = 1mol/L	

Method:	Basic Experimental Parameters:	
	Titrant delivery rate (mL/min.)	6
	No. of endothermic endpoints	1
	Data smoothing factor	40
	ERC (2 nd derivative) value	-12
	Procedure:	
	Prepare a titration vessel with an deionized water and tare on a balance approximately 10mL of industrial phose and weigh. Record the mass of the sam until thoroughly mixed. If necessary, of to near room temperature. Titrate with 1 mol/L solution to a single thermometry	pproximately 30mL e. Slowly pipette in sphoric acid sample mple. Carefully swirl cool in a water bath standard $c(BaCl_2) =$ ric endpoint.
	It is good practice to remove the titrat thermometric sensor immediately aft finished, rinsing the titration assemb water.	ion vessel from the er each titration is bly thoroughly with
	Blank determination. Weigh accurately 6, 8 and 10g of a typical sample of in acid into titration vessels containing 3 titrate as above. Plot sample mas endpoint volume (y-axis), and determin regression analysis. This volume is the "blank" value, to be subtracted from the in subsequent determinations. R correlation coefficient for this determinations.	approximately 2, 4, ndustrial phosphoric 30mL DI water, and as (x-axis) against be the y-intercept by a systematic error or the endpoint volume ecord the linear ation.

A Metrohm

Results (Example):	Its (Example): Analysis of sulfate in works phosphoric acid:		icid:
		As SO4 ²⁻ % w/w	1
	Mean (n=10):	2.52±0.01	
	Blank determination	Blank = 0.1522mL	$R^2 = 0.9993$

Calculation:	% SO ₄ = <u>(EP mL- blank mL)*c(BaCl₂)*96.0626*0.1)</u>	
	Sample size, g	

