


# Thermo. Titr. Application Note No. H-013

**Title:** Determination of Moisture Content of Ultra-fine Solids

**Scope:** Determination of water in moist particulate solids such as cobalt oxyhydroxide.

**Principle:** A direct titration of the moisture content of ultra-fine or of high porosity solids may yield inconsistent and low results due to the inability of the titrant to rapidly penetrate moisture-containing interstices. In such cases, a back-titration may be recommended. An excess of 2,2-dimethoxypropane (DMP) is added to a suspension of the moisture-containing solid in a suitable polar solvent such as propan-2-ol to which a small amount of a suitable acid catalyst has been added. After thorough mixing of the suspension with the excess DMP, this excess is back-titrated with a standard solution of water in propan-2-ol. Water reacts endothermically with DMP to form acetone and methanol.



**Reagents:** 2mol/L 2,2-dimethoxypropane (DMP) in dry (HPLC grade) cyclohexane [DMP data: FW = 104.15, d = 0.847, purity = 98%]. Dilute 250mL DMP to 1000mL with cyclohexane in a volumetric flask  
Methane sulfonic acid  
Propan-2-ol (HPLC grade)  
2 mol/L H<sub>2</sub>O in propan-2-ol

**Method:** Basic Experimental Parameters:  
Data rate (per second) 10  
Titrant delivery rate (mL/min.) 2  
No. of exothermic endpoints 1  
Data smoothing factor 75

<b>Method (Continued):</b>	Weigh rapidly and accurately approximately 0.5 - 1g moist sample (depending on moisture content) into a dry titration vessel. Deliver by bulb pipette 25.00mL dry propan-2-ol, plus 250µL methane sulfonic acid into the vessel.  Start the titration. The pre-dose part of the back-titration commences. 5.00 mL DMP reagent is dispensed into the titration vessel. A 30-second delay is programmed prior to the back-titration with standard H <sub>2</sub> O in propan-2-ol so that all moisture is extracted from the sample and reacted with the DMP. The titration proceeds to an endothermic endpoint.  The DMP reagent is standardized against the standard water solution by dispensing volumes of DMP ranging from 1-5 mL, and titrating them with 25mL propan-2-ol and 250µL methane sulfonic acid.
----------------------------	---

<b>Results:</b>	Analysis of very fine moist cobalt oxyhydroxide			
	<b>Sample Mass, g</b>	<b>DMP Pre-Dose, mL</b>	<b>Std. H<sub>2</sub>O Back Titre, mL</b>	<b>% H<sub>2</sub>O</b>
0.5102	5.000	1.921	18.9	
0.5067	5.000	1.977	18.7	
0.5030	5.000	2.002	18.6	
0.5180	5.000	1.865	19.0	
0.4906	5.000	2.035	18.9	
0.4934	5.000	2.032	18.8	
0.5306	5.000	1.829	18.8	
<b>Average</b>				<b>18.8</b>
<b>Standard Deviation</b>				<b>0.1</b>


***Calculation***

Data Item	Symbol	Test Data
DMP pre-dose vol, mL	P <sub>D</sub>	5.000
DMP molarity	M <sub>D</sub>	1.9610
Std. H <sub>2</sub> O soln. molarity	M <sub>W</sub>	2.0001
Std. H <sub>2</sub> O back-titre, mL	T <sub>W</sub>	2.032
Titration blank, mL	T <sub>O</sub>	0.3001
Sample mass, g	S	0.4934
Formula weight H <sub>2</sub> O	F <sub>W</sub>	18.015

$$\% \text{ } H_2\text{O} = \frac{(((P_D \times M_D) - (M_W(T_W + T_O))) \times F_W \times 100)}{(S \times 1000)}$$

$$\% \text{ } H_2\text{O} = \frac{(((5.000 \times 1.9610) - (2.0001 \times (2.032 + 0.3001))) \times 18.015 \times 100)}{(0.4934 \times 1000)}$$

$$= 18.8\%$$

**Thermometric Titration Plot:****Legend:**

Red = solution temperature curve

Black = second derivative curve