

## Thermo. Titr. Application Note No. H-018

Title:	Determination of Chloride in Drilling Fluids	
Scope:	Determination of chloride in oil well drilling fluids	
Principle:	Acidification of the sample with 2mol/L HNO <sub>3</sub> and titration with 0.1M AgNO <sub>3</sub> . Some emulsified samples with heavy o loadings may require addition of a hydrocarbon solven such as mineral turpentine to achieve adequate phase separation.	
Reagents:	0.1M AgNO <sub>3</sub> solution	
	Concentrated HNO <sub>3</sub>	
	Mineral turpentine or kerosene	
Method:	Basic Experimental Parameters:	
	Data rate (per second) 10	
	Titrant delivery rate (mL/min.)	
	No. of exothermic endpoints 1	
	Data smoothing factor 40	
	Procedure: It has been found that with some drilling fluids, the matrix is very oily and immiscible with the acidified fluid in the beaker. Low results with poor reproducibility can be obtained.	
	Weigh accurately approximately 1g of thoroughly homogenized drilling fluid into a clean, dry titration vessel. Add a spin bar, 10mL of mineral turpentine (obtainable from a hardware store) 10mL of deionized water, and 5mL 2 mol/L HNO <sub>3</sub> . Titrate to a single exothermic endpoint with 0.1M AgNO <sub>3</sub> .	

Results:			
	Sample	Chloride as %CaCl <sub>2</sub> w/v	
	А	3.05, 2.99, 2.85	
	В	3.89, 3.85, 3.84	
	С	2.13, 2.13, 2.11	
	D	1.87, 1.89, 1.94	

1



CaCl<sub>2</sub> 
$$g / L = \frac{((Titre - blank) \times FW \ CaCl_2 \times M \ AgNO_3 \times 100)}{(sample \ mass, \ g \times 1000)}$$

$$CaCl_2 \ g / L = \frac{((Titre - blank) \times FW \ CaCl_2 \times M \ AgNO_3 \times 100)}{(sample \ mass, \ g \times 1000)}$$

