Thermo. Titr. Application Note No. H-045

 Title:
 Standardization of EDTA titrant by magnesium

| Scope: | Standardization of tetrasodium EDTA titrant for use in the |
|--------|--|
| | determination of magnesium. |

| Principles: | Tetrasodium EDTA (Na ₄ EDTA) is the preferred reagent for the thermometric complexometric titration of metals, due to its much higher solubility than the normally used dibasic salt Na ₂ H ₂ EDTA. |
|-------------|---|
| | The thermometric titration of magnesium with EDTA is carried out in an ammonia/ammonium chloride buffer (~pH 10) environment. The endpoint is marked by a slight upswing in temperature, as the reaction between EDTA and magnesium is endothermic. |

| Reagents: | Titrant. 1 mol/L tetrasodium EDTA |
|-----------|--|
| | Buffer: NH_3/NH_4CI solution, pH 10. Dissolve 70g NH_4CI in 688mL conc. NH_3 soln. and make to 1000mL with D.I. water. |
| | Standard magnesium solution. Sufficient A.R. magnesium metal (ribbon) to make 500mL of a 0.2 mol/L solution is scaped clean of an oxide to a uniform shiny surface and degreased. After weighing, it is transferred to a wide mouth erlenmeyer flask, and covered with 100mL D.I. water. 25mL conc. HCl is added slowly through a funnel in the mouth of the flask, maintaining a steady effervescence until all the metal is dissolved. The solution is cooled, and quantitatively transferred to a 500mL volumetric flask, making to volume with D.I. water. |

| Method: | Basic Experimental Parameters: | | | |
|---------|---|---|--|--|
| | Titrant delivery rate (mL/min.) | 2 | | |
| | No. of endpoints | 1 | | |
| | Data smoothing factor | 50 | | |
| | Stirring speed (802 stirrer) | 6 | | |
| | Delay before start (secs.) | 15 | | |
| | Buffer pre-dose (from Dosino), mL | 5 | | |
| | Pipette aliquots of Mg standard soluvessel. Allowing for the addition of 5m make up the difference in the volum water. Titrate to change in gradient curve. | aliquots of Mg standard solution into a titration Allowing for the addition of 5mL of buffer solution, p the difference in the volume to 30mL with DI Fitrate to change in gradient in the temperature | | |

| Results (example): | | | | | |
|---|-------------|------------------------|--------------|--|--|
| For this exercise, ~99% Mg (Lab. Reagent grade) was used. | Aliquot, mL | mmole Mg ²⁺ | Titre, mL | | |
| | 25 | 4.9893 | 4.941, 4.947 | | |
| | 20 | 3.9914 | 3.965, 3.965 | | |
| | 15 | 2.9936 | 2.982, 2.981 | | |
| | 10 | 1.9957 | 2.001, 2.001 | | |
| | 5 | 0.9979 | 1.004, 1.007 | | |
| | | | | | |



