

Thermo. Titr. Application Note No. H-047

| | |
|---------------|---|
| Title: | Determination of Nickel by EDTA Back-Titration |
|---------------|---|

| | |
|---------------|--|
| Scope: | Determination of nickel in refinery and plating solutions. When other metals capable of being complexed by EDTA are present, these will interfere and enhance the result for nickel. |
|---------------|--|

| | |
|--------------------|---|
| Principles: | <p>Tetrasodium EDTA (Na_4EDTA) is the preferred reagent for the thermometric complexometric titration of metals, due to its much higher solubility than the normally used dibasic salt $\text{Na}_2\text{H}_2\text{EDTA}$.</p> <p>Nickel reacts too slowly with EDTA to be considered for a direct thermometric titration. In this determination, an excess of EDTA is added, sufficient time allowed to complete the reaction, and the excess EDTA back-titrated with standard Cu solution.</p> <p>The thermometric titration of nickel (and associated metals) with EDTA is carried out in an ammonia/ammonium chloride buffer (~pH 10) environment. The endpoint is marked by a slight upswing in temperature, caused by the formation of the Cu ammine complex, whose reaction is slightly more exothermic than that of the Cu-EDTA complex</p> |
|--------------------|---|

| | |
|------------------|--|
| Reagents: | <p>Titrant 1. 1 mol/L tetrasodium EDTA (standardized) – Dosino #1</p> <p>Titrant 2. 1 mol/L Cu^{2+} solution, prepared from A.R. $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ and standardized against standard Na_4EDTA solution – Dosino #2</p> <p>Buffer: $\text{NH}_3/\text{NH}_4\text{Cl}$ solution, pH 10. Dissolve 70g NH_4Cl in 688mL conc. NH_3 soln. and make to 1000mL with D.I. water – Dosino #3</p> |
|------------------|--|

| | | |
|---|------------------------------------|----|
| Method: | Basic Experimental Parameters: | |
| | Titrant delivery rate (mL/min.) | 2 |
| | No. of exothermic endpoints | 1 |
| | Data smoothing factor | 60 |
| | Stirring speed (802 stirrer) | 6 |
| | Delay before start (secs.) | 60 |
| | Buffer pre-dose (mL) | 5 |
| | Na ₄ EDTA pre-dose (mL) | 5 |
| For this exercise, a 0.15 mol/L Ni ²⁺ solution was prepared from A.R. NiSO ₄ ·6H ₂ O. A 20mL aliquot is pipetted into a titration vessel, and the titration program commenced. Titrate to a change in gradient of the temperature curve. | | |

| | |
|---------------------|--|
| Calculation: | |
| | $\%Ni = \frac{(((EDTA \text{ mL} - (Cu^{2+} \text{ mL} \times MCu^{2+} / MEDTA)) \times MEDTA \times FW Ni \times 100)}{(Sample \text{ mass, g} \times 1000)}$ |

| | |
|---|------------------------|
| Results (example): | |
| Sample: A.R. NiSO ₄ ·6H ₂ O | %Ni = 99.99±0.14, n=11 |

Thermometric Titration Plot:

Legend:

Red = solution temperature curve

Black = second derivative curve

