Thermo. Titr. Application Note No. H-048

Title:	Standardization of disodium	
	dimethylglyoximate for the determination of	
	nickel	

Scope:	Standardization	of	disodium	dimethylgloximate	by
	thermometric titrat	tion v	vith standard	d Ni(II) solution.	

Principle:	Standardization	of	disodium	dimethylglyoximate
	(Na ₂ DMG) solution	by titr	ation with sta	ndard Ni solution to
	an exothermic en	dpoint	in buffered	ammonia solution.
	Two moles of DMG	Freact	with one mol	e of Ni(II). Acidic Ni
	solutions possibly	contail	ning Al(III) an	id Fe(III) should be
	complexed with tar	trate p	rior to basifica	ation.

Reagents:	Titrant: 0.5 mol/L disodium dimethylglyoximate. Dissolve 153.6 g disodium dimethylglyoximate (99% pure, FW = 304.21) in deionized water and make to 1000mL in a volumetric flask.
	NH_3/NH_4CI buffer: Dissolve 17.5 g A.R. NH_4CI in 172 mL A.R. conc. NH_3 soln. Make to 250 mL with deionized water.
	Kna tartrate solution, 450g/L.
	Standard Ni solution, 0.1 mol/L. In the exercise described here, a reference sample of nickel oxide was kindly supplied by a nickel refinery. The stated purity was 78.4 \pm 0.24% Ni, determined by gravimetric analysis by the dimethylgloxime method. Sufficient NiO to make 500mL of 0.1 mol/L solution was weighed into a 150mL beaker provided with a magnetic stirrer. 20mL concentrated HNO ₃ and 60mL HClO ₄ was added and heated with stirring until near boiling, continuing the heating until all the oxide was dissolved (a little colloidal silica remained). The solution was cooled and made to 500mL with DI water in a volumetric flask.

Method:	hod: Basic Experimental Parameters:			
	Titrant delivery rate (mL/min.)	2		
	No. of exothermic endpoints	1		
	Data smoothing factor	60		
	Stirrer speed	12*		
	Procedure: The sodium of dispensed from a 10mL Dosino. was also dispensed from a Dosir buffer solution using the pre- Titrotherm software.	limethylglyoximate was The standard Ni solution no, as was the NH ₃ /NH ₄ Cl e-dose function of the		
	If a Dosino is used to dispense the is only necessary to add 30mL of tartrate solution to the titration determination. The tartrate solut interference from Al(III) and Fe(containing Ni, and is added in the solution to "matrix match". The volume of 40-45mL prior to the normal volume together with reduces the detrimental effect nickel dimethylglyoximate precipi	If a Dosino is used to dispense the standard Ni solution, it is only necessary to add 30mL of DI water and 5mL KNa tartrate solution to the titration vessel prior to the determination. The tartrate solution is used to prevent interference from Al(III) and Fe(III) impurities in samples containing Ni, and is added in the case of high purity Ni solution to "matrix match". The aim is to have a solution volume of 40-45mL prior to the titration. This higher than normal volume together with the high stirrer speed* reduces the detrimental effect that the highly viscous nickel dimethylglyoximate precipitate has on the endpoint.		
	Equipment stained by the precipitate may be cleaned by s strong mineral acid; eg, 1 standardization exercise, the cleaned between each titration, a DI water.	Ni dimethylglyoximate soaking in a solution of a mol/L HCI. During the titration assembly was and rinsed thoroughly with		

Results:	Titration of aliquots of standard 0.10036 mol/L Ni(II) solution, dispensed by Dosino			
	mL Ni(II) soln.	mmole Ni	mL Na ₂ DMG soln.	
	2.50	0.2503	1.050	
	3.75	0.3755	1.526	
	5.00	0.5007	2.028	
	6.25	0.6259	2.533	
	7.50	0.7510	3.043	
	8.75	0.8762	3.541	



