

Thermo. Titr. Application Note No. H-058

Title:	Determination of Sodium as Chloride in Ketchup and Sauces			
	·			
Scope:	Determination of sodium as chloride in ketchups, sauces and similar food products.			
Principle:	Titration with standard AgNO ₃ solution to a single exothermic endpoint. Sodium content is determined assuming it is all present as NaCl.			
Reagents:	Titrant: 0.1 mol/L AgNO ₃ solution. pH modifier: 5 mol/L HNO ₃ solution			
Method:	Basic Experimental Parameters:			
	Titrant delivery rate (mL/min.) 2			
	No. of exothermic endpoints 1			
	Data smoothing factor (DSF) 50			
	Stirring speed (802 stirrer) 7			
	Stir time before titration starts (sec.) 25			
	Sample preparation:			
	Depending on estimated sodium content of the sample, weigh 0.5 – 2g of sample into a clean, dry titration vessel. Add 30mL D.I. water and 1mL 5 mol/L HNO ₃ . Adjust sample mass to obtain a titre of ~3 – 5mL 0.1 mol/L AgNO ₃ .			
	Note: Since food labelling regulations demand the result be expressed in <i>mg Na/100mL</i> (rather than <i>mg Na/100g</i>), it is necessary to separately determine the density of the product.			



Examples:				
	Brand	Claimed total Na mg/100mL (on label)	Chloride content, expressed as Na g/100mL	
	Master Foods tomato sauce	1141	833±2.1 (n=5)	
	Master Foods tomato sauce, reduced salt	598	536±1.2 (n=5)	
	Fountain tomato sauce	835	689±1.3 (n=5)	
	Fountain barbecue sauce	550	435±1.3 (n=5)	
	Rosella tomato sauce	1250	997±0.3 (n=7)	
	Coles tomato sauce	910	811±0.4 (n=7)	

$\label{eq:calculations:} \begin{aligned} \text{Na g/100g} &= \frac{((\text{Titre, mL} - \text{blank, mL}) \times \text{AgNO}_3 \ \text{mol/L} \times 22.98977 \times 100}{\text{sample mass, g}} \\ \text{Na g/100mL} &= \frac{\text{Na g/100g}}{\text{density}} \end{aligned}$

