Thermo. Titr. Application Note No. H-081

Title: Determination of Phosphoric and Nitric Acid in Nitrophos Liquors

Scope:	Determination	of	phosphoric	and	nitric	acids	in	liquors
	from the Nitrop	hos	s fertilizer ma	nufa	cturing	g proce	SS.	

Principle:	Titration in a standardized solution of 2 mol/L NaOH in a
	sodium chloride solution to sharpen the third endpoint.

Reagents:	<i>Titrant:</i> 2 mol/L NaOH, standardized against potassium
	hydrogen phthalate
	Conditioning solution: 180g/L NaCl

Method:	Basic Experimental Parameters:	
	Titrant delivery rate (mL/min.)	5
	Delay before titration commences (se	c) 5
	No. of exothermic endpoints	3
	Data smoothing factor (DSF)	93
	Stirring speed (802 stirrer)	15
	<i>Titration:</i> Weigh accurately ap (approximately 1g) of Nitrophos liqu titration vessel. Add 30mL 180g/L Nato stand for approximately 30 minut order that any contained hexafluor hydrolyzed.	proximately 0.7mL for into a clean, dry Cl solution, and allow es before titrating in rosilicic acid is fully
	Blank determinations: Prepare titration weighing masses of Nitrophos lique range from approximately 1.5 to 0.6g two graphs. In the first graph, plot se against titration volume to the find Determine the y-intercept from regree will be the blank (B1) for the nitric at the second graph, plot sample mass difference between the third and sec EP2. Once again, determine the regression analysis. This will be the phosphoric acid determination	n solutions as above, for encompassing a . Titrate, and prepare sample mass (x-axis) irst endpoint, EP1. ession analysis. This icid determination. In s (x-axis) against the cond endpoints, EP3- te y-intercept from e blank (B2) for the

Examples:	Nitrophos process liquor		
	H ₃ PO ₄ % w/w	24.1±0.25 (n=11)	
	HNO ₃ % w/w	7.8±0.17 (n=11)	

Calculations:	
	$H_{2}PO_{4} \% = \frac{((EP3 - EP2 - B2) \times NaOH mol/L \times 97.9769 \times 100)}{(EP3 - EP2 - B2) \times NaOH mol/L \times 97.9769 \times 100)}$
	(sample mass, $g \times 1000$)
	HNO $\%$ = (((EP1 - B1) - (EP3 - EP2 - B2)) × NaOH mol/L × 62.99564 × 100)
	(sample mass, $g \times 1000$)

