

# Thermo. Titr. Application Note No. H-114

**Title:** Determination of sulfuric acid, nitric acid, and hydrofluoric acid in etch solutions

**Scope:** Determination of sulfuric acid, nitric acid, and hydrofluoric acid etch solutions.

**Principle:** Two separate titration sequences are required to analyze the mixture:

- titration of the HF content with  $\text{Al}(\text{NO}_3)_3$  (the “elpasolite” reaction)
- titration of the  $\text{H}_2\text{SO}_4$  with  $\text{BaCl}_2$  followed by titration with NaOH to determine the “total acids” content

The HF,  $\text{H}_2\text{SO}_4$ , and “total acids” contents are converted to a  $\text{HNO}_3$  equivalent, with the  $\text{HNO}_3$  content found by subtracting the HF and  $\text{H}_2\text{SO}_4$  from the “total acids” content.

**Reagents:**

*Titrant #1:* 0.5 mol/L  $\text{Al}(\text{NO}_3)_3$  in D.I. water  
*Titrant #2:* 1 mol/L  $\text{BaCl}_2$  in D.I. water  
*Titrant #3:* 2 mol/L NaOH  
*Standard #1:* Anhydrous A.R. NaF, 99%  
*Standard #2:* Anhydrous A.R.  $\text{Na}_2\text{SO}_4$ , 99.9%  
*Standard #3:* A.R. potassium hydrogen phthalate, 99.95%  
*pH 4.5 buffer solution.* Weigh 130.9 g anhydrous potassium acetate, and 54.7 g anhydrous acetate and dissolve in 600 mL deionized water. Add 115 mL glacial acetic acid, and make to 1000 mL with deionized water.

**Method:**

*Basic equipment list (for automated titration):*  
859 Titrotherm interface  
814 Sample Processor with 24 place rack, 75 mL PP tubes  
802 Magnetic Stirrer with propeller mixer  
800 Dosinos, 1 x 5 mL, 3 x 10 mL burettes  
6.9011.040 Thermoprobe (fluoride resistant)

*Basic experimental parameters:*

*HF titration*

|                                       |     |
|---------------------------------------|-----|
| Sample size (g):                      | ~3  |
| Titrant delivery rate (mL/min)        | 4   |
| ERC (exothermic)                      | -10 |
| Data smoothing factor                 | 53  |
| Stirring speed (802 Magnetic Stirrer) | 12  |

|  |  |      |    |
|--|--|------|----|
| <i>H<sub>2</sub>SO<sub>4</sub> titration</i> |  |      |    |
| Titrant delivery rate (mL/min)               |  | 4    |    |
| ERC (exothermic)                             |  | -25  |    |
| Data smoothing factor                        |  | 50   |    |
| Stirring speed (802 Magnetic Stirrer)        |  |      | 14 |
| <br>   |  |      |    |
| <i>“total acids” titration</i>               |  |      |    |
| Titrant delivery rate (mL/min)               |  | 4    |    |
| ERC (exothermic)                             |  | -150 |    |
| Data smoothing factor                        |  | 50   |    |
| Stirring speed (802 Magnetic Stirrer)        |  | 15   |    |

*Methods:*

*NOTES:*

- (a) *In each of the following methods, the HF titration must precede the H<sub>2</sub>SO<sub>4</sub>-“total acids” titration*
- (b) *Due to the concentrated nature of the solution, use of an “air pipette” to dispense accurate aliquots is not recommended. In the case of an accurate auto-diluter not being available, it is permissible to weigh in the samples, correcting the sample mass to a volume equivalent via a previously determined density of the sample.*

*1. HF titration*

When installing the method for the first time, create a CV to hold the result for HF g / L. Dispense accurately 2 mL of sample solution into a plastic titration vessel, add 5 mL D.I. water, and add 30 mL acetate buffer, either automatically via a Dosino or by auto dilution. Titrate with the HF titration method.

*2. H<sub>2</sub>SO<sub>4</sub>/total acids titrations*

Dispense accurately 1 mL of sample solution into a plastic titration vessel, and add 30 mL D.I. water. Titrate with the H<sub>2</sub>SO<sub>4</sub> “total acids” method.

*3. Determination of blank for HF titration*

When installing the method for the first time, prepare a CV for the HF titration blank. Prepare at least 5 plastic titration vessels, with dispensing volumes ranging from 0.5 to 2 mL of etch solution. Add 5 mL D.I. water and 30 mL acetate buffer. Titrate with the HF titration blank method.

*4. Determination of blanks for H<sub>2</sub>SO<sub>4</sub> and “total acids” titrations*

When installing the methods for the first time, create CVs for the H<sub>2</sub>SO<sub>4</sub> and “total acids” blanks. Prepare at least 5 plastic titration vessels, dispensing volumes of sample

ranging from 0.2 to 1.0 mL. Titrate with the H<sub>2</sub>SO<sub>4</sub> "total acids" titration blank method.

*5. Standardization of 0.5 mol/L Al(NO<sub>3</sub>)<sub>3</sub> titrant*

Prepare at least 5 plastic titration vessels, weighing in accurately approximately 0.25 to 0.6 g freshly dried anhydrous NaF in roughly even increments. Add 5 mL acetate buffer and 25 mL water. Titrate with the Al(NO<sub>3</sub>)<sub>3</sub> standardization method.

*6. Standardization of 1 mol/L BaCl<sub>2</sub> titrant*

Prepare at least 5 plastic titration vessels, weighing in accurately approximately 0.14 to 0.85 g freshly dried anhydrous A.R. Na<sub>2</sub>SO<sub>4</sub> in roughly even increments. Add 5 mL 5 mol/L HNO<sub>3</sub> and 25 mL D.I. water. Titrate using the BaCl<sub>2</sub> standardization method.

*7. Standardization of 2 mol/L NaOH titrant*

Prepare at least 5 plastic titration vessels, weigh in accurately approximately 0.4 to 2.4 g freshly dried A.R. potassium hydrogen phthalate in roughly even increments. Add 30 mL D.I. water, and titrate using the NaOH standardization method.

**Results:**

*Analysis of a synthetic H<sub>2</sub>SO<sub>4</sub>-HNO<sub>3</sub>-HF mixture*

|                      | <i>H<sub>2</sub>SO<sub>4</sub></i> | <i>HF</i> | <i>HNO<sub>3</sub></i> |
|----------------------|------------------------------------|-----------|------------------------|
| <i>Mean, g/L</i>     | 430.4                              | 37.1      | 287.9                  |
| <i>Std. dev. g/L</i> | 0.5                                | 0.1       | 0.7                    |
| <i>% RSD</i>         | 0.11                               | 0.22      | 0.24                   |

**Calculations:**

*HF titration:*

$$\text{HF g/L} = ((\text{EP mL} - \text{HF blank, mL}) * \text{Al(NO}_3)_3 \text{ mol/L} * \text{FW HF} * 6) / \text{sample volume, mL}$$

*Conversion to HNO<sub>3</sub> g/L equivalent:*

$$\text{HNO}_3 \text{ g/L equivalent} = \text{HF g/L} * \text{FW HNO}_3 / \text{FW HF}$$

*H<sub>2</sub>SO<sub>4</sub> "total acids" titrations*

$$\text{H}_2\text{SO}_4 \text{ g/L} = ((\text{EP mL} - \text{H}_2\text{SO}_4 \text{ blank, mL}) * \text{BaCl}_2 \text{ mol/L} * \text{FW H}_2\text{SO}_4) / \text{sample volume, mL}$$

*Conversion to HNO<sub>3</sub> g/L equivalent:*

**Calculations**

(continued)

$$\text{HNO}_3 \text{ g/L equivalent} = \text{H}_2\text{SO}_4 \text{ g/L} \cdot 2 \cdot \text{FW HNO}_3 / \text{FW H}_2\text{SO}_4$$

$$\text{“total acids” (as HNO}_3\text{) g/L} = ((\text{EP mL} - \text{“total acids” blank, mL}) \cdot \text{NaOH mol/L} \cdot \text{FW HNO}_3 / \text{volume, mL})$$

**Titration Plots:**




Fig. 1. Example of HF with Al titration plot




Fig. 2. Example of H<sub>2</sub>SO<sub>4</sub> with BaCl<sub>2</sub> titration plot

**Titration plots  
(continued)**




Fig. 3. Example of „total acids“ with NaOH titration plot