


Determination of ferrous ion in acidic solutions with ceric solution






This Application Note looks at the determination of ferrous ion in acidic solutions from approximately 0.25 g/L by thermometric titration with ceric titrant. The exothermic oxidation reaction shows a sharp endpoint that is detected using the Thermoprobe as a sensitive temperature sensor.

Method description

Principle

Ceric ion in strongly acidic solution is a suitable titrant for the direct determination of ferrous ion.



Thus, one mole Ce^{4+} corresponds to one mole Fe^{2+} . The reaction is exothermic. Although solution temperature rises are not as large as with $\text{Cr}_2\text{O}_7^{2-}$ or MnO_4^{-} , endpoints are sharp. Ce^{4+} solutions have the advantage over MnO_4^{-} solutions in that they are indefinitely stable, and only require standardization on preparation. If desired, perchloric acid may be substituted for sulfuric acid to enhance the redox potential.

E° volts in 1 mol/L acid solutions:

in H_2SO_4 : 1.44

in HClO_4 : 1.70

Samples

«Test solution»

Sample preparation

No sample preparation required

Configuration

Basic equipment list for automated titration

814 USB Sample Processor	2.814.0030
859 Titrotherm	2.859.0010
Sample rack 24 × 75 mL	6.2041.340
Thermoprobe, fluoride resistant	6.9011.040
Sample beaker 75 mL	6.1459.400
802 Rod Stirrer	2.802.0010
Stirring propeller 104 mm	6.1909.020
1 × 800 Dosino	2.800.0010
1 × Dosing unit 10 mL	6.3032.210
tiamo™	6.6056.222

Solutions

Titrant	$c(\text{Ce}^{4+}) = 0.25 \text{ mol/L}$ ammonium ceric nitrate, $(\text{NH}_4)_2\text{Ce}(\text{NO}_3)_6$ in 1 mol/L H_2SO_4 or HClO_4 solution
Acid	10% (w/v) H_2SO_4 solution, or 20% (w/v) HClO_4 solution
Standard solution	$c(\text{Fe}^{2+}) = 0.3 \text{ mol/L}$ ferrous ammonium sulfate (FAS), $(\text{NH}_4)_2\text{SO}_4 \cdot \text{FeSO}_4 \cdot 6\text{H}_2\text{O}$, prepared in 0.1 mol/L H_2SO_4 solution. Prepare freshly for each batch of Ce^{4+} solution to be standardized. Calculate amount of FAS to be weighed, taking into account the stated purity of the reagent.
Note: FAS is not a primary standard, but is suitable for standardizing Ce^{4+} solutions used for routine analysis. It is advised that Ce^{4+} solutions should be prepared in large quantities (10 to 25 L) to avoid the necessity of frequent standardizations.	
«Test solution»	A test solution of approx. ~1 g/L Fe^{2+} was prepared from FAS to simulate a customer solution containing a low Fe^{2+} concentration. This was prepared by weighing approximately 1.77 g FAS into a 250 mL volumetric flask containing 10 mL 10% (w/v) H_2SO_4 and making to volume with dist. water.

Analysis

Basic method (for approx. 1 g/L Fe^{2+} solution)

Pipette a 20 mL aliquot of test solution into a titration vessel, add 10 mL 10% (w/v) H_2SO_4 , and titrate with $c(\text{Ce}^{4+}) = 0.25 \text{ mol/L}$ solution to an exothermic endpoint.

Standardization

Using a specially-constructed **tiamo™** program, aliquots from 0.5 to 5.0 mL were dispensed in 0.5 mL increments from a Dosino into titration vessels containing 10 mL of 10% (w/v) H_2SO_4 and volumes of dist. water to give a total volume of approx. 30 mL at the start of titration. The **tiamo™** program performs the regression analysis and automatically calculates the molarity, titration blank, and correlation coefficient of the regression line.

Method description

Parameters

Titrant dose rate (mL/min)	4
ERC EP1 (exothermic)	-15
Data smoothing («filter factor»)	45
Stirring speed (802 Rod Stirrer)	14
Evaluation start (mL)	0.5
Damping until (mL)	0.5

Calculations

$$\text{g/L Fe}^{2+} = ((\text{EP1} - \text{blank}) \times \text{C01} \times \text{C02}) / \text{C00}$$

EP1 = titrant consumption in mL

C00 = sample weight in mL


C01 = concentration of cer titrant in mol/L

C02 = molecular weight of Fe (55.847 g/mol)

Results

«test solution» Fe^{2+}	1.03 ± 0.005 g/L (theoretical = 1.006 g/L)
Standardization:	c(Ce^{4+}) = 0.2389 mol/L blank = 0.029 mL coefficient of determination, $R^2 = 1.0000$

Titration plot



Titration plot, Fe^{2+} with Ce^{4+}