

# ISE Application Note No. I - 3

<b>Title:</b>	<b>Sodium content of electrolyte powder</b>
<b>Summary:</b>	Determination of sodium in electrolyte powder (for pharmaceutical use) by direct potentiometry using the Sodium ISE
<b>Sample:</b>	Pharmaceutical product (electrolyte powder)
<b>Sample Preparation:</b>	Exactly weigh about 0.5 g sample into a 100 mL volumetric flask, dissolve in dist. water, fill up to the mark and mix.
<b>Instruments and Accessories:</b>	692 pH/Ion Meter, 725 Dosimat, 728 Magnetic Stirrer, 6.0501.100 Na ISE, 6.0726.100 Ag/AgCl reference electrode (inner electrolyte c(HCl) = 0.1 mol/L, outer electrolyte TISAB), printer
<b>Reagents:</b>	<b>Sodium standard:</b> 1000 mg/L Na <sup>+</sup> , prepared from NaCl (2.542 g/L NaCl) <b>TISAB:</b> Mix 75 mL triethanolamine with 500 mL dist. water. Set the pH value to 8.0 by adding HCL (c(HCl) = 2 mol/L). Fill up to 1000 mL with dist. water.
<b>Analysis:</b>	<b>Calibration:</b> Fill 25 mL of TISAB into the measuring vessel, then start the automatic calibration with a 5 point calibration curve. Values should be between 0.001 and 0.01 mol/L Na <sup>+</sup> . <b>Sample:</b> Add 25 mL TISAB and 1.0 mL sample solution into the measuring vessel, then start the automatic determination by pressing the «print» button of the 692 pH/Ion Meter.
<b>Remarks:</b>	After each determination the electrode has to be rinsed with dist. water. The electrode should then be «dried» carefully with a soft paper tissue. <b>Result:</b> AVG(5) = 20.22 +/- 0.2 % Na <sup>+</sup>