

ISE Application Note No. I - 8

Title:	Sulphide content of waste water
Summary:	Determination of sulphide in waste water by direct potentiometry with the Silversulphide ISE
Sample:	Different waste water samples
Sample Preparation:	(see remarks)
Instruments and Accessories:	692 pH/Ion Meter, 725 Dosimat, 728 Magnetic Stirrer, 6.0502.180 Ag ₂ S ISE, 6.0726.100 Ag/AgCl reference electrode (c(KCl) = 3 mol/L), printer
Reagents:	Sulphide standard: 5000 mg/L S ²⁻ , prepared from Na ₂ S * 9 H ₂ O (37.46 g/L). Titre has to be determined by potentiometric titration with AgNO ₃ . Antioxidant ISA: Contains 2 mol/L NaOH and 0.2 mol/L each, ascorbic acid and Na ₂ EDTA.
Analysis:	Pipette 50 mL ISA into a beaker and deaerate it with nitrogen. Add 10.0 mL sample. Start the automatic determination by three standard additions and a preselected delta U of 12 mV.
Remarks:	To avoid loss of volatile sulphides in neutral or acidic solutions add 10 mL c(NaOH) = 10 mol/L to the sample bottle (1 L) prior to sampling. Replace the outer electrolyte of the reference electrode daily. For low sample concentrations use lower standard concentrations. Results: sample a) AVG(3) = 67.2 +/- 0.2 mg/L S ²⁻ sample b) AVG(3) = 73.3 +/- 1.1 mg/L S ²⁻ sample c) AVG(3) = 225 +/- 0.5 mg/L S ²⁻