KF Application Note No. K-4

Title: Water in lyophilisates (e.g. vaccines in vials)

Summary:	The water content of lyophilisates contained in vials is determined			
according to Karl Fischer. Conditioned solvent (methano				
	into the vial to dissolve the sample and extract the water (ultrasonic			
	bath). Afterwards the contents of the vial are transferred to the titra-			
	tion vessel to carry out the automatic determination.			

Sample:	Three different lyophilised vaccines in vials		
Sample			
Preparation:	none		

Instruments and Accessories:	701 KF Titrino or 720 KFS Titrino, 703 Titration Stand, printer, ultra- sonic bath		
Analysis:	Pour 25 mL methanol into the titration vessel and start conditioning. When the solution has been conditioned take out 10 mL with a sy- ringe and inject 5 mL of this solution through the septum into the vial. Dissolve the sample by putting the vial into the ultrasonic bath. Transfer the contents of the vial to the titration vessel using a syringe and start the automatic titration.		
	Reagents:Solvent:methanol (dry)Titrant:Hydranal Composite 2 (Riedel-de Haën)		
Results:	With the described procedure the total water content of a vial is de-		

Results: With the described procedure the total water content of a vial is determined. To calculate the water content of the product alone one has to know the sample size.

Settings:	701 KF Titrino		
	>titration parameter	S	
	extr.time	0 s	
	stop crit.:	drift	
	stop drift	20 uL/min	
	>preselections		
	conditioning:	on	
	report:	full	