KF Application Note No. K-5

Water in ink	
The water content of in	k is determined according to Karl Fischer.
Ink (main component: t	ripropylene glycol monoethyl ether)
none	
1	
701 KF Titrino or 720	KFS Titrino, 703 Titration Stand, printer
, , ,	the exact sample mass is determined by dif-
e e ,	n a syringe into the titration vessel containing
	nt (methanol). The solvent has to be changed
	JH.
Reagents:	
-	rv)
(omposite 5 (Riedel-de Haën)
,	
AVG(3) = 0.326 +/- 0.0	0006 % water
701 KF Titrino	
>titration parameters	
extr.time	0 s
•	drift
•	25 uL/min
•	a n
•	on
	on full
	The water content of in Ink (main component: to none 701 KF Titrino or 720 I Inject ca. 3 g sample (ference weighing) with the conditioned solver after each determination Reagents: Solvent: methanol (co Titrant: Hydranal Co AVG(3) = 0.326 +/- 0.00 701 KF Titrino >titration parameters