KF Application Note No. K-15

Title:	Water in sweet liquorice	
Summary:	Fischer. To dissolve the	weet liquorice is determined according to Karl he sample a mixture of methanol and formam- and a high-frequency mixer as stirring device.
Sample:	Sweet liquorice	
Sample Preparation:	Cut up the sample into small pieces.	
Instruments and		
Accessories:	701 KF Titrino or 720 KFS Titrino, 703 Titration Stand, printer, Poly- tron PT 1200 Disintegrator	
Analysis:	Pour 50 mL solvent mixture into the titration vessel and condition it. Wait for a steady drift below 10 ug/min. Add ca. 200 mg sample us- ing a weighing spoon (the exact sample mass is determined by dif- ference weighing). Start the PT 1200 Disintegrator during the extrac- tion time. Carry out at least a triplicate determination.	
	Reagents:	
	•	xture: volume ratio methanol : formamide =
	Titrant: Hydranal C	omposite 5 (Riedel-de Haën)
Decultor		
Results:	AVG(8) = 15.45 +/- 0.3	38 % Waler
Settings:	701 KF Titrino	
	>titration parameters	
	extr.time	120 s
	stop crit.:	drift
	stop drift	20 uL/min
	>preselections	
	conditioning:	on
	req.smpl size:	ON full
	report:	full