KF Application Note No. K-26

Title:	Water in used lubricating oil	
Summary:	The water content of used lubricating oil is determined according to Karl Fischer by coulometric titration. To prevent unwanted side reactions special KF reagents are used.	
Sample:	Used lubricating oil	
Sample		
Preparation:	none	
	,	
Instruments and Accessories:	737 KF Coulometer, cell with diaphragm; 728 Magnetic Stirrer, printer	
Analysis:	Start the 737 KF Coulometer to condition the reagent in the meas- uring cell. When the drift is steady and below 10 ug/min inject ca. 0.4 g sample using a syringe. Exchange the reagent in the cathode compartment as soon as ist colour changes from brown to yellow.	
	Reagents: Hydranal Coulomat CK and Hydranal Coulomat AK (Riedel-de Haën)	
Results:	AVG(9) = 631 +/- 23 ppm water	

Settings:	737 KF Could	737 KF Coulometer		
	smpl.req. d.start extr. stop drift: delay time report:	on 16 ug/min 10 s auto 3 s full		