

KF Application Note No. K-29

Title: Water in spices (marjoram, nutmeg, pepper)

Summary: The water content of spices is determined according to Karl Fischer.

To release the water from the cells a high-frequency mixer has to be

used.

Sample: Dried marjoram, ground nutmeg and pepper

Sample

Pepper corns and nutmegs have to be cooled prior to grinding. Preparation:

Instruments and

701 KF Titrino, 720 KFS Titrino or 758 KFD Titrino, 703 Titration Accessories:

Stand, printer, Polytron PT 1200 Disintegrator

Pour ca. 50 mL dry methanol into the the titration vessel and condi-Analysis:

> tion it. Add ca. 0.2 ... 0.4 g sample using a weighing spoon and switch on the PT 1200 Disintegrator for 120 s (extraction time), then carry out the water determination. Exchange the methanol after two

determinations.

Reagents:

Solvent: methanol (dry)

Titrant: Hydranal Composite 5 (Riedel-de Haën)

Results: Marjoram: AVG(5) = 8.46 + -0.06 % water

> AVG(5) = 6.07 + -0.06 % water Nutmeg: White pepper: AVG(5) = 8.36 + -0.07 % water

Settings: 701 KF Titrino

>titration parameters

extr.time 120 s stop crit.: drift

stop drift 20 uL/min

>preselections

conditioning: on req.smpl size: on full report: