

KF Application Note No. K- 44

Title: Water in panthenol

Summary: The water content in panthenol is determined according to Karl Fischer

Sample: Panthenol

Sample Preparation: Weigh about 1 g of the viscous sample into a dry Erlenmeyer flask and add about 15 g dry methanol. Close the flask with a septum stopper and mix.

Instruments and Accessories: 701 or 720 or 758 Titrino, 703 Titration Stand, Printer

Analysis: Add about 20 mL Working Medium and about 1 mL Composite K to the titration vessel and condition it. Add about 2 g prepared sample mixture using a syringe and titrate the water content. Do the same for a blank with the methanol.

Reagents:

Solvent: Hydranal Working Medium K (Riedel de Haen)

Titrant: Hydranal Composite 5K (Riedel de Haen)

Results: AVG(3) = 3.97 +/- 0.10 % water

Settings:

720 KFS-Titrino		> control parameters	
> titration parameters		EP at U	250 mV
titr. direction:	„-“	dynamics	100 mV
pause 1	0 s	stop crit:	drift
Start V:	off	stop drift	20 uL/min
extr. Time	0 s		
I(pol)	50 uA		