
IC Application Note No. N-2

Title:	Methylarsonic acid and Dimethylarsinic acid	
Summary:	Determination of methyarsonic acid and dimethylarsinic acid using anion chromatography with direct conductivity detection	
Sample:	Standard solution	
Sample Preparation:	-	
Column:	6.1006.100 Metrosep Anion Dual 2	
Eluent:	1 mmol/L phosphate buffer, pH = 5.3	
Flow:	0.8 mL/min	
Injection Volume:	100 µL	



Results:	Methylarsonic acid mg/L	Dimethylarsinic acid mg/L
	10	10