

IC Application Note No. N-49

Title:	Borate and silicate in washing powder
--------	---------------------------------------

Summary: Determination of borate and silicate in washing powder using anion chromatography with direct conductivity detection.

Sample: Washing powder

Sample Preparation: Dissolution and subsequent dilution 1: 10000, injection

through an IC-RP cartridge


Column: 6.1005.000 Hamilton PRP-X 100

Eluent: 3.2 mmol/L sodium hydroxide,

0.5 mmol/L sodium carbonate

Flow: 1.0 mL/min

Injection Volume: 20 µL



Results:	Borate g/kg	Silicate g/kg
	59.7	34.7