IC Application Note No. N-51

Title:	concentra	•	ate and sulfa sion solutior otides		
Summary:	Determination of acetate, chloride, citrate and sulfate in a concentrate of an infusion solution using anion chromatography with direct conductivity detection. Non-suppressed IC is used to avoid interferences by the amino acids.				
Sample:		Concentrate of an infusion solution containing: 75 g/L amino acids, 30 g/L dipeptides, 1g/L citric acid			
Sample Preparation:	Dilution 1 : 5	0 with eluent			
Column:	6.1009.000 Super-Sep Anion				
Eluent:	2.5 mmol/L p	2.5 mmol/L phthalic acid, 2% acetonitrile $pH = 4.2$ (TRIS)			
Flow:	1.5 mL/min	1.5 mL/min			
Injection Volume:	20 µL				
uS/cm 152- 151- 150- 149- 148- 147- 0 1 2	Chloride Chloride	-Sulfate		System Peak	
Results:	Acetate g/L	Chloride g/L	Citrate g/L	Sulfate g/L	
Sample	2.710	1.600	1.095	n.q.	