IC Application Note No. N-56

Title:	Silicate besides fluoride in tap water		
Summary:	Determination of silicate besides fluoride in tap water using anion chromatography with direct conductivity		

detection.

Sample: Tap water

Dilution 1: 1000 with ultrapure water Sample Preparation:


6.1005.100 Phenomenex Starlon A 300 (2 in a row) Column:

6.1005.020 PRP-X 100 Guard for preconcentration

Eluent: 1.0 mmol/L sodium hydroxide, 15 mmol/L fructose

Flow: 1.5 mL/min

Injection Volume: 10 mL preconcentrated



Results:	Silicate mg/L	Fluoride mg/L
	6.42	n.q. [*]

Comment: The eluent was purged with helium during operation. The fluoride concentration of Herisau tap water is approx. 40 µg/L.