
# Identification of 45 aromatic and medicinal plants used in cosmetic and pharmaceutical industry



This Application Note shows that Vis-NIR spectroscopy can identify 45 different aromatic and medicinal plants (e.g. Organicum majoricum and Tilia cordata). The characteristic spectrum of each plant, represented in a spectral library, can be used in raw material identification. In contrast to microscopic identification — which takes several hours and requires trained chemists — the presented Vis-NIR method is extremely fast and can be performed by unskilled operators.



# Method description

### Introduction

The global market potential for medical and aromatic plants is on increase in the recent past. Medical plants are used as raw materials for production of pharmaceuticals, cosmetics, and food supplements. Therefore, they are not only the basic compounds of traditional medicines, but also modern medicines are derived from medicinal plants. Aromatic plants provide the raw material for the production of flavors in e.g., personal care products. The quality control of medical plants is usually done by microscopic identification or wet chemistry of plant characteristics. Vis-NIR spectroscopy can be used to replace the time-consuming conventional identification method.

### Configuration

The following equipment was used (Tab. 1).

Tab. 1:

Equipment	Metrohm order code
NIRS DS2500 Analyzer	2.922.0010
Large sample cup	6.7402.050
Vision 4.03 Software	6.6069.102

## **Experimental**

A NIRS DS2500 analyzer (reflection measurement), was used to collect the spectral data of multiple samples of 46 different medical and aromatic plants (**Fig. 1**). Sample heterogeneity was handled using a rotating sample cup, which allows averaging of 8 single scans. In sum, the natural occurring variation is represented in 546 representative spectra.



Fig. 1: A NIRS DS2500 analyzer was used to record the representative spectra.



In Vision (Metrohm chemometrical software), using the algorithm of residual variance with a threshold of 0.9, a library for the following substances (Tab.2), was developed. The entire wavelength range, covering Visand NIR range (400 – 2500 nm), was used and the data were pretreated using a  $2^{\rm nd}$  derivative combined with a Standard Normal Variate (SNV). Cross-validation was applied to verify the performance of the derived library.

Tab. 2

Aloysia citriodora
Angelica archangelica
Angelica montana
Arctium lappa
Calendula officinalis
Camelia sinensis
Cinnamomum verum
Crataegus monogyna
Crocus sativus
Curuma longa
Elataria cardamomum
Eleutherococcus
senticosus
Equisetum arvense
Euphrasia stricta
Fraxinus excelsior
Galanga officinalis
Gentiana lutea
Glycyrhiza glabra
Harpagophytum
procumbens
Hibiscus sabadarifa
Hypericum perforatum
Illiciums verum

Lavendula angustifolia
Lavendula intermedia
Matricaria recutita
Meliisa officinalis
Mentha piperita
Organicum majorana
Origanicum hirtum
Origanum majoricum
Plantago lanceolata
Rosa damascene
Rosmarinus officinalis
Sambucus nigra
Schinus mole
Schinus mole
Schinus mole Thymus vulgaris
Schinus mole Thymus vulgaris Tilia cordata
Schinus mole Thymus vulgaris Tilia cordata Tilia platyphyllos
Schinus mole Thymus vulgaris Tilia cordata Tilia platyphyllos Tilia tomentosa
Schinus mole Thymus vulgaris Tilia cordata Tilia platyphyllos Tilia tomentosa Valeriana officinalis

### Results

The characteristic spectrum of 45 different aromatic and medicinal plants, represented in a spectral library can be used in raw material identification. **Fig. 2** shows the example spectrum of 4.




Fig. 2: Characteristic spectra of Organicum majoricum, Tilia cordata, Illicum verum and Harpagophytum procumbens