


# Sorbitol and sucrose in soap applying pulsed amperometric detection



Polyols and sugar components in soaps increase the stability and size of bubbles and foam. A typical soap ingredient is glycerol. But sugar alcohols and sugars are used as well. The determination of sorbitol and sucrose in a transparent soap is achieved on a Metrosep Carb 2 - 150/4.0 column with subsequent pulsed amperometric detection (PAD).

## Results

Compound	Concentration [g/kg]	RSD [% , n = 3]	Recovery [%]
Sorbitol	77.7	0.0	100
Sucrose	269.4	0.6	102

## Sample

Transparent soap

## Sample preparation

0.1 g of soap is dissolved in 1000 mL of ultrapure water.  
Injection after filtration (0.2 µm).

## Columns

Metrosep Carb 2 - 150/4.0	6.1090.420
Metrosep Carb 2 Guard/4.0	6.1090.500

## Solutions

<u>Eluent</u>	100 mmol/L sodium hydroxide 10 mmol/L sodium acetate
---------------	---

## Parameters

Flow rate	0.5 mL/min
Injection volume	5 µL
P <sub>max</sub>	20 MPa
Recording time	20 min
Column temperature	32 °C

## PAD Parameters


Cell	Wall-Jet cell
Working electrode	Gold
Reference electrode	Palladium
Spacer	50 µm
Measuring potential	0.05 V
Measuring duration	100 ms
Cycle duration	550 ms
Measuring range	200 µA
Temperature	35 °C
Mode	PAD

## Analysis

Pulsed amperometric detection

## Instrumentation

940 Professional IC Vario ONE	2.940.1100
IC Amperometric Detector	2.850.9110
919 IC Autosampler plus	2.919.0020
IC equipment Wall-Jet cell: Carb (Au, Pd)	6.5337.010



[www.metrohm.com](http://www.metrohm.com)

 **Metrohm**