
# Stability Application Note R–012 Oxidation stability of different solid foodstuffs



The 892 Professional Rancimat can be used for determining the oxidation stability of fats contained in solid, non-meltable foodstuffs such as peanuts, almonds, potato chips (crackers), cookies, french fries. The crushed samples can be used directly without extraction of the fat.

Resu	lts

Sample	Induction time
Peanuts	3.1 h
Almonds	7.6 h
Potato chips (crackers)	28.5 h
Butter cookies	20.6 h
French fries	6.0 h



# Method description

## Sample

Various types of solid fat-containing foodstuffs

#### Sample preparation

Samples were crushed into small pieces. Powdered samples can be used, too.

French fries were cut into small pieces with a knife.

### Instrument

892 Professional Rancimat



#### Parameters

Temperature	120 °C
Air flow	20 L/h
Evaluation	Induction time
Evaluation sensitivity	1.0
Stop criteria	End points
Sample amount	0.5 g
Measuring solution	50 mL deionized water

www.metrohm.com

