


# IC Application Note No. S-142

**Title:** Silicate, sulfate and phosphate in a clay extract

**Summary:** Determination of silicate, sulfate and phosphate in clay extract using anion chromatography with conductivity detection before and after chemical suppression. Applying a step gradient and valve switching to work with or without chemical suppression.


**Sample:** Clay extract  
**Sample preparation:** Dilution 1 : 200, injection after filtration

**Column:** 6.1005.100 Phenomenex Starlon A 300  
**Eluent A:** 1.0 mmol/L sodium hydroxide, 0.1 mmol/L sodium carbonate  
**Eluent B:** 1.0 mmol/L sodium hydroxide, 1.0 mmol/L sodium carbonate  
**Suppressor:** MSM (50 mmol/L H<sub>2</sub>SO<sub>4</sub>)  
**Flow:** 2.0 mL/min  
**Injection Volume:** 100 µL



<b>Results:</b>	Silicate as SiO <sub>2</sub> mg/kg	Sulfate mg/kg	Phosphate mg/kg
		1970	498

Switching from non-suppressed to suppressed mode



Valve B:

Time 0.0 ... 4.5 min INJECT position, i.e. non-suppressed mode for silicate  
polarity = -

Time 4.5 ... 30 min FILL position, i.e. suppressed mode for sulfate and phosphate  
polarity = +