

IC Application Note No. S-183

Title: Advanced inline dialysis setup for Ion Chromatography

Summary: Determination of acetate, chloride, nitrate, phosphate and sulfate in mayonnaise using anion chromatography with conductivity detection after chemical suppression and advanced dialysis as inline sample preparation.

Sample: Mayonnaise

Sample Preparation: 0.5 g dissolved in 100 g ultrapure water, injection after inline dialysis


Column: 6.1006.510 Metrosep A Supp 5 – 100

Eluent: 3.2 mmol/L sodium carbonate
1.0 mmol/L sodium hydrogen carbonate

Suppressor: Metrohm Suppressor Module (MSM, 50 mmol/L H₂SO₄)


Flow: 0.7 mL/min

Injection Volume: 20 µL after inline dialysis



Results:	Acetate g/L	Chloride g/L	Nitrate	Phosphate mg/L	Sulfate mg/L
	9.012	7.186	n.q.	68.3	44.6

Advanced Dialysis Setup:



Advanced inline dialysis compared to the former setup

Old System:

Table with RSD of 40 injections of 5 ppm anions:

Ion	Mean	StdDev	RSD%	Recov.
Fluoride	5.105	0.084	1.6	102.1
Chloride	5.196	0.084	1.6	103.9
Nitrite	5.260	0.077	1.5	105.2
Bromide	5.203	0.067	1.3	104.1
Nitrate	5.139	0.076	1.5	102.8
Phosphate	4.737	0.168	3.5	94.8
Sulfate	4.868	0.114	2.3	97.4

Table with carry over of the old and new system:

Ion	old system in %	new system in %
Fluoride	0.22	0.12
Chloride	0.25	0.14
Nitrite	0.12	0.08
Nitrate	0.21	0.01
Phosphate	0.25	0
Sulfate	0.34	0.25

New System:

Table with RSD of 40 injections of 5 ppm anions:

Ion	Mean	StdDev	RSD%	Recov.
Fluoride	4.900	0.040	0.8	98.0
Chloride	4.936	0.041	0.8	98.7
Nitrite	5.029	0.044	0.9	100.6
Bromide	4.968	0.047	1.0	99.4
Nitrate	4.920	0.028	0.6	98.4
Phosphate	4.847	0.066	1.4	96.9
Sulfate	4.847	0.057	1.2	96.9

Comment:

After an injection of a 5 ppm standard, a blank (ultrapure water) was injected. This was done 5 times to calculate the carry over.

Old system:

Stopped flow achieved by means of a 6 port valve.