

IC Application Note No. S-232

Title: Chloride and sulfate in electrolytes used in sensors for transcutaneous CO₂ measurement

Summary: Determination of chloride and sulfate in an electrolyte used in sensors for transcutaneous CO₂ measurement applying anion chromatography with conductivity detection after chemical suppression.

Sample: Electrolyte used in sensors for transcutaneous CO₂ measurement

Sample Preparation: Dilution 1 : 500 with ultrapure water, followed by inline dialysis


Column: 6.1006.510 Metrosep A Supp 5 – 100

Eluent: 3.2 mmol/L sodium carbonate
1.0 mmol/L sodium hydrogencarbonate

Suppressor: Metrohm Suppressor Module (MSM, 50 mmol/L sulfuric acid)

Flow: 0.7 mL/min

Injection Volume: 20 µL



Results:

Chloride
mg/L / mol/L

5761 / 0.16

Sulfate
mg/L

21.5