

# IC Application Note No. S-253

**Title:** Molybdate in 2.5% NaCl applying Inline Matrix Elimination by sample re-injection

**Summary:** Determination of molybdate in 2.5% NaCl using anion chromatography with conductivity detection after chemical suppression and Inline Matrix Elimination by molybdate preconcentration after the first separation and subsequent re-injection.

**Sample:** 1 mg/L molybdate in 2.5% NaCl

**Sample Preparation:** Direct injection

**Column:** 6.1005.350 Metrosep A Supp 1 HS


**Preconc. Column:** 6.1006.300 Metrosep A PCC 1

**Eluent:** 3.0 mmol/L sodium carbonate

**Suppressor:** Metrohm Suppressor Module (MSM, 200 mmol/L H<sub>2</sub>SO<sub>4</sub>)

**Flow:** 1.5 mL/min

**Injection Volume:** 20 µL




**Results:**

Molybdate  
mg/L

1

System setup:



Step 1: Injection of the sample on «Valve 2»

Step 2: Molybdate collection on the preconcentration column on «Valve 1»

Step 3: Re-injection of the collected molybdate with «Valve 1»

Molybdate is collected onto concentrating column

