IC Application Note S-304

Variable Inline Preconcentration including Matrix Elimination for trace anion analysis (MiPCT-ME)



Metrohm Inline Preconcentration Technique with Matrix Elimination (MiPCT-ME) is a powerful tool that combines preconcentration, matrix elimination, and multilevel calibration. The latter only requires a single multi-ion standard solution. The 800 Dosino takes over all liquid handling tasks. The shown system setup allows sample analysis from 0.1 μ g/L up to 1.0 mg/L.

Results

	Conc. [ng/L]		Conc. [ng/L]
Fluoride	100	Nitrate	100
Chloride	100	Phosphate	100
Nitrite	100	Sulfate	100
Bromide	100		



Sample

Standard solution

Sample preparation

Inline Preconcentration with Matrix Elimination (MiPCT-ME)

Columns

Metrosep A Supp 5 - 150/4.0	6.1006.520
Metrosep A Supp 4/5 Guard/4.0	6.1006.500
Metrosep A PCC 1 HC/4.0	6.1006.310

Solutions

<u>Eluent</u> (inline eluer preparation)	nt 3.2 mmol/L sodium carbonate 1.0 mmol/L sodium hydrogen carbonate
Suppressor regenerant	100 mmol/L sulfuric acid
Rinsing solution	Ultrapure water

Analysis

Conductivity after sequential suppression

Parameters

Flow rate	0.7 mL/min
Injection volume	404000 µL
P _{max}	15 MPa
Recording time	18 min
Column temperature	35 °C

Instrumentation

850 Professional IC Anion – MCS	2.850.2030
IC Conductivity Detector	2.850.9010
858 Professional Sample Processor	2.858.0010
2 x 800 Dosino (liquid handling)	2.800.0010
849 Level Control for Inline Eluent Preparation	2.849.1030

Calibration MiPCT-ME

Calibration range	Factor of 100
Standard solution:	
All anions	10 μg/L
1. Level	40 μL = 0.1 μg/L
2. Level	80 μL = 0.2 μg/L
3. Level	200 μL = 0.5 μg/L
4. Level	400 μ L = 1.0 μ g/L
5. Level	800 μL = 2.0 μg/L
6. Level	2000 μL = 5.0 μg/L
7. Level	4000 μ L = 10 μ g/L



www.metrohm.com



Flow chart



Procedure: The Dosino aspirates the required sample volume into the transfer tubing. Then the sample volume is loaded on to the preconcentration column, which subsequently is rinsed with excess ultrapure water.

Due to the fact that the sample volume can be widely adjusted, the calibration from $0.1...10 \mu g/L$ converts into a calibrated measuring range from 0.1 up to 1000 $\mu g/L$.

Injecting 40 μL of a 1000 $\mu g/L$ sample equals injecting 4000 μL of a 10 $\mu g/L$ solution (the highest calibration level).

www.metrohm.com

