


# Trace anions on Metrosep A Supp 16 - 250/2.0 after Inline Preconcentration and Matrix Elimination (MiPCT-ME)



Intelligent Preconcentration Technique with Matrix Elimination (MiPCT-ME) is used for trace determination of the seven standard anions. On a microbore Metrosep A Supp 16 - 250/2.0 column, the analysis is concluded within 33 minutes with high recovery rates. The detection limits are in the low ng/L range (calculated by MagIC Net) using a preconcentration volume of 2000 µL.

## Results

Anion 0.5 µg/L	Recovery [%]	LOD [µg/L]	Anion 0.5 µg/L	Recovery [%]	LOD [µg/L]
1 Fluoride	97	0.001	5 Nitrate	99	0.005
2 Chloride	98	0.002	6 Sulfate	99	0.003
3 Nitrite	98	0.003	7 Phosphate	97	0.012
4 Bromide	99	0.005			

### Sample

0.5 µg/L QC standard

### Sample preparation

Intelligent Preconcentration Technique with Matrix Elimination (MiPCT-ME)

### Columns

Metrosep A Supp 16 - 250/2.0	6.1031.230
Metrosep A Supp 16 Guard/2.0	6.1031.600
Metrosep A PCC 2/4.0	6.1006.330

### Solutions

Eluent	7.5 mmol/L sodium carbonate 0.75 mmol/L sodium hydroxide
Suppressor regenerant	100 mmol/L <u>sulfuric acid</u>
Rinsing solution	STREAM

### Analysis

Conductivity detection after sequential suppression

### Parameters


Flow rate	0.25 mL/min
Injection volume (MiPCT-ME)	2000 µL
P <sub>max</sub>	16 MPa
Recording time	33 min
Column temperature	45 °C

### Instrumentation

940 Professional IC Vario ONE/SeS/PP	2.940.1500
IC Conductivity Detector	2.850.9010
858 Professional Sample Processor	2.858.0010
2 x 800 Dosino	2.800.0010
MSM Rotor A	6.2832.000
Adaptor sleeve Suppressor Vario	6.2842.020
IC equipment: MiPCT-ME	6.5330.160



### Schematic setup



[www.metrohm.com](http://www.metrohm.com)

 **Metrohm**