


Metrohm IC Driver for OpenLab CDS: Anions with Dose-in Gradient and Dosino Regeneration

Working with Dosinos in OpenLab: Determination of anions applying a Dose-in Gradient, Dosino Regeneration, and full control by OpenLab CDS.



Anion chromatogram of eight anions while applying a Dose-in gradient, acquired by OpenLab CDS.

OpenLab CDS is the newest generation of chromatography data systems from Agilent. The Metrohm IC Driver 1.0 for OpenLab CDS implements Metrohm ion chromatographs in OpenLab CDS for full control and data acquisition. This application shows the use of a gradient (Dose-in Gradient) as well as Dosino Regeneration in OpenLab CDS. Fluoride, chloride, nitrite, bromide, nitrate, sulfate, phosphate, and iodide in a standard solution are separated and determined.

Results

Anion	Result [mg/L]	Anion	Result [mg/L]
1 Fluoride	1.0	5 Nitrate	5.0
2 Chloride	5.0	6 Sulfate	5.0
3 Nitrite	1.0	7 Phosphate	5.0
4 Bromide	1.0	8 Iodide	5.0

Sample

Standard solution

Sample preparation

Direct injection

Columns

Metrosep A Supp 17 - 150/4.0	6.01032.420
Metrosep A Supp 17 Guard/4.0	6.01032.500

Solutions

Eluent A	5.0 mmol/L sodium carbonate 0.2 mmol/L sodium hydrogen carbonate
Eluent B (Dosino)	50 mmol/L sodium carbonate 2.0 mmol/L sodium hydrogen carbonate
Regenerant	100 mmol/L phosphoric acid
Rinsing	STREAM

Instrumentation

930 Compact IV Flex SeS/Deg	2.930.1460
IC Conductivity Detector	2.850.9010
858 Professional Sample Processor	2.858.0020
800 Dosino (gradient)	2.800.0010
800 Dosino (MSM regeneration)	2.800.0010
MSM Rotor A	6.2832.000
Adapter sleeve for Suppressor Vario	6.2842.020
IC equipment: Dosino regeneration	6.5330.190
IC equipment: Additional eluent	6.5330.090
EMPTY sensor for 2 L bottles	6.2769.110
Metrohm IC Driver 1.0 for OpenLab CDS	6.6080.100
OpenLab CDS 2.4 (Agilent)	

Analysis


Conductivity detection after sequential suppression

Parameters

Flow rate	0.7 mL/min
Injection volume	20 µL
P _{max}	18 MPa
Column temperature	Ambient
Recording time	19 min



Dose-in Gradient



www.metrohm.com

 **Metrohm**