


IC Application Note S-386

Fluoride in sodium fluoride and phosphoric acid gel for pharmaceutical use

USP monograph modernization: ion chromatography applying a hydroxide eluent on a Metrosep A Supp 16 - 250/4.0 column (L91)



Chromatogram of the System Suitability test for assay (concentrations; sodium fluoride 2 $\mu\text{g/mL}$, sodium acetate 1 $\mu\text{g/mL}$). Insert: chromatogram of the assay (corresponding to 1 $\mu\text{g/mL}$ sodium fluoride).

Sodium fluoride and phosphoric acid gel for pharmaceutical use need to comply with USP requirements. The actual monograph (USP 42) uses two different methods for the identification and the assay. Ion chromatography allows the measurement of these two parameters within a single determination. In the course of the USP monograph modernization, this ion chromatographic approach makes this type of analysis even easier.

Results

Anion	Sample ID	Result [%]	USP limit [%]
1 Fluoride	Assay [%]	99.8	90–110
2 Acetate	Impurity	n.q.	-
3 Chloride	Impurity	n.q.	-

n.q. = not quantified. For further results, see next page.

Sample

Dental gel containing sodium fluoride and phosphoric acid.

Sample preparation

Dilution with ultrapure water. Final concentration 1 µg/mL.

Columns

Metrosep A Supp 16 - 250/4.0	6.1031.430
Metrosep A Supp 16 Guard/4.0	6.1031.500

IC Solutions

Eluent	15 mmol/L potassium hydroxide
Regenerant Dosino	500 mmol/L sulfuric acid
Rinsing	Ultrapure water

Instrumentation

930 Compact IC Flex Oven/ChS/PP/Deg	2.930.2360
IC Conductivity Detector	2.850.9010
858 Professional Sample Processor	2.858.0020
800 Dosino for Dosino Regeneration	2.800.0010
M5M-HC Rotor A	6.2842.000
IC equipment: Dosino Regeneration	6.5330.190

Parameters IC

Flow rate	1.0 mL/min
Injection volume (MiPT)	20 µL
P _{max}	20 MPa
Column temperature	40 °C
Recording time	27 min

Analysis

Conductivity detection after chemical suppression

System suitability requirements for assay

USP Parameter	Result	USP required	Remarks
Resolution F ⁻ /acetate	NLT 1.5	5.4	Pass
Tailing factor F ⁻	NMT 2.0	1.1	Pass
RSD F ⁻ [%; n=5)	NMT 2%	0.24	Pass

NLT = not less than, NMT = not more than

