

Ti Application Note No. T- 31

Title: **Na₂O (free base) and SiO₂ (silicate) in water glass**

Summary: Determination of Na₂O and SiO₂ in water glass by potentiometric titration with HCl using the Sb electrode.

Sample: Water glass (solution of Na silicate)

Sample Preparation: none

Instruments and Accessories: 702, 716 or 736 Titrino or 726 Titroprocessor, 6.0303.000 Sb electrode and 6.0712.100 reference electrode (consists of 6.0724.140 internal system + 6.1227.000 electrolyte vessel)

Analysis: Weigh exactly ca. 1 g sample into a plastic beaker and dissolve in 50 mL dist. water. Titrate with c(HCl) = 1 mol/L and stop after the first equivalence point (free base). Add 5 mL potassium fluoride solution (400 g/L), wait for 2 min and titrate again with c(HCl) = 1 mol/L (silicate).

Calculation: $\% \text{Na}_2\text{O} = \text{C31} = \text{EP1} * \text{C01} * \text{C02} * \text{C04} / (\text{C00} * \text{C03})$

excess HCl = C30 = C41 - EP1

EP1 = titrant consumption in mL for the first titration

C00 = ca. 1 (sample weight in g)

C01 = 1 (concentration of the titrant in mol/L)

C02 = 62 (M(Na₂O) in g/mol)

C03 = 20 (special factor)

C04 = titre of the titrant

C41 = titrant consumption in mL at the end of the first titration

Calculation: % SiO₂ = RS1 = (EP1 + C30) * C01 * C02 * C03 * C04 / C00

Mass ratio = RS1 / C31

Solids = RS1 + C31

EP1 = titrant consumption in mL for the second titration

C00 = ca. 1 (sample weight in g)

C01 = 1 (concentration of the titrant in mol/L)

C02 = 60.1 (M(SiO₂) in g/mol)

C03 = 0.025 (special factor)

C04 = titre of the titrant

C30 = excess HCl in mL from the first titration (Na₂O)

C31 = result % Na₂O (first titration)

Remarks: Polish the surface of the Sb electrode from time to time with wetted Alox powder.

Results:

AVG(6) = 9.53 +/- 0.02 % Na₂O

AVG(6) = 30.00 +/- 0.05 % SiO₂

AVG(6) = 39.53 +/- 0.05 % solids

AVG(6) = 3.15 +/- 0.01 (mass ratio)