

Ti Application Note No. T- 36

Title:	Benzydamine hydrochloride in disinfectant solution
---------------	---

Summary:	Determination of benzydamine hydrochloride {1-benzyl-3-[3-(dimethylamino)-propoxy]-1H-indazole hydrochloride} in disinfectant solution by potentiometric titration with sodium tetraphenylborate using the NIO Surfactant Electrode.
-----------------	--

Sample:	Disinfectant solution
----------------	-----------------------

Sample Preparation:	none
----------------------------	------

Instruments and Accessories:	702, 716, 736 or 751 Titrino or 726 Titroprocessor, 6.0507.010 NIO Surfactant Electrode and 6.0733.100 Ag/AgCl reference electrode
-------------------------------------	--

Analysis:	Pipette 10.0 mL sample solution into a beaker, add ca. 90 mL dist. water and 4 ... 5 drops glacial acetic acid and titrate with sodium tetraphenylborate $c(\text{STPB}) = 0.01 \text{ mol/L}$.
------------------	--

Calculation:	$\text{mg/mL benzydamine hydrochloride} = \text{EP1} * \text{C01} / \text{C00}$ EP1 = titrant consumption in mL C00 = 10.0 (sample size in mL) C01 = 3.4586 (benzydamine hydrochloride equivalent in mg/mL; 1 mL $c(\text{STPB}) = 0.01 \text{ mol/L}$ corresponds to 3.4586 mg benzydamine hydrochloride)
---------------------	---

Remarks:	Result: AVG(3) = 1.583 +/- 0.002 mg/mL benzydamine hydrochloride
-----------------	---