

# Ti Application Note No. T- 43

**Title:** Bromine index in low-level standards

**Summary:** Determination of the bromine index in low-level standards by bi-voltametric titration with bromide/bromate using a double Pt electrode.

**Sample:** Two different standard solutions

**Sample Preparation:** none

**Instruments and Accessories:** 702, 716, 736 or 751 Titrino or 726 Titroprocessor, 6.0308.100 double Pt wire electrode, titration vessel with thermostatic jacket, cryostat

**Analysis:** Put 100 mL solvent mixture\* into the titration vessel and add the sample (the sample size depends on the bromine index). Wait with stirring until the temperature is between 0 ... 5 °C, then titrate with  $c(\frac{1}{6} \text{KBrO}_3) = 0.02 \text{ mol/L}$  (also containing 2.04 g/L KBr) using the MET Ipol mode with a polarisation current of 1 uA.

\*) Solvent mixture:  
714 mL glacial acetic acid + 134 mL 1,1,1-trichloroethane + 134 mL methanol + 18 mL  $w(\text{H}_2\text{SO}_4) = 20 \%$

**Calculation:** Bromine index (= mg bromine / 100 g)  
 $= (\text{EP1} - \text{C01}) * \text{C02} * \text{C03} * \text{C04} / \text{C00}$

EP1 = titrant consumption for the sample in mL  
C00 = sample mass in g  
C01 = titrant consumption for the blank (solvent) in mL  
C02 = 0.02 (equivalent concentration of the titrant in mol/L)  
C03 = 7990 (100 \* M(Br) in g/mol)  
C04 = titre of the titrant (determined with cyclohexene)

**Remarks:** Results:  
Sample A: 85.7 mg bromine / 100 g  
Sample B: 9.5 mg bromine / 100 g