

Ti Application Note No. T- 45

Title: Acetate, chloride and phosphate in infusion solution

Summary: Determination of acetate, chloride and phosphate in infusion solution by potentiometric titration with sodium hydroxide after conversion of the anions to the corresponding acids.

Sample: Infusion solution

Sample Preparation: Pass 50 mL solution through a strongly acidic cation exchanger resin (column) into a 100 mL volumetric flask. Rinse with dist. water, fill the flask to the mark and mix the contents.

Instruments and Accessories: 702, 716, 736 or 751 Titrino or 726 Titroprocessor, 6.0233.100 combined pH glass electrode

Analysis: Pipette 50.0 mL of the prepared sample solution (corresponds to 25 mL of the original sample) into a beaker and titrate with $c(\text{NaOH}) = 0.2 \text{ mol/L}$.

The first equivalence point of the titration curve corresponds to the sum of HCl and the first dissociation step of H_3PO_4 , the second EP corresponds to acetic acid and the third EP to the second dissociation step of H_3PO_4 .

Calculation: 1 mL $c(\text{NaOH}) = 0.2 \text{ mol/L}$ corresponds to
11.689 mg NaCl
19.629 mg CH_3COOK
27.217 mg KH_2PO_4

$$\text{RS1} = \text{EP3} - \text{EP2}$$

$$\text{RS2: g/L KH}_2\text{PO}_4 = \text{RS1} * \text{C01} / \text{C00}$$

$$\text{RS3: g/L NaCl} = (\text{EP1} - \text{RS1}) * \text{C02} / \text{C00}$$

$$\text{RS4: g/L CH}_3\text{COOK} = (\text{EP2} - \text{EP1}) * \text{C03} / \text{C00}$$

Calculation:	EP1 = titrant consumption to reach the first EP
	EP2 = titrant consumption to reach the second EP
	EP3 = titrant consumption to reach the third EP
	C00 = 25.0 (mL of original sample contained in the sample volume used for the titration)
	C01 = 27.217 (KH ₂ PO ₄ equivalent in mg/mL)
	C02 = 11.689 (NaCl equivalent in mg/mL)
	C03 = 19.629 (CH ₃ COOK equivalent in mg/mL)

Remarks:	<ul style="list-style-type: none">– If the sample also contains sulphate, this is detected together with HCl and the first dissociation step of H₃PO₄ at the first EP of the titration curve. In this case chloride has to be determined separately by potentiometric titration with AgNO₃.– Gluconate can not be differentiated from acetate by titration.
-----------------	---