

Ti Application Note No. T- 56

Title: Lauryl ether sulphate by photometric/turbidimetric titration

Summary: Determination of lauryl ether sulphate (LAES) by potentiometric/turbidimetric titration with TEGO trant A100 using the 610 nm Spectrode.

Sample: Lauryl ether sulphate standards

Sample Preparation: none

Instruments and Accessories: 702, 716, 736, 751 or 785 Titrino, 722 Propeller Rod Stirrer, 6.5501.01X Spectrode 610 nm, Metrodata TiNet 2

Analysis: Pipette the sample into a beaker, add 10 mL buffer solution pH = 3 and make up to about 60 mL with dist. water, then titrate with c(TEGO trant A100) = 0.004 mol/L using the MET mode.

Calculation: 1 mL c(TEGO trant A100) = 0.004 mol/L corresponds to 0.004 mmol LAES

$$\text{mmol/L LAES} = \text{EPx} * \text{C01} * \text{C02} / \text{C00}$$

EPx = titrant consumption in mL
EP1: evaluation Titrino
EP2: intersection point TiNet

C00 = sample volume in mL
C01 = 0.004
C02 = 1000 (conversion factor in mL/L)

Remarks: This method cannot be used for formulations containing non-ionic surfactants.

Results:	Date	11.08.1998	Time	11:09:46
	User	JS		
	Method	LAES with Spectrode		
	Id1	610nm MET		
	EPC	30		
	SmplSize	4 ml		
	Endpoints:			
	MET U.EP1	9.395 ml	322 mV	
	MET U.Intersection1	9.496 ml	0 mV	
	Results:			
LAES1	9.395 mmol/l			
LAES2	9.496 mmol/l			