

Ti Application Note No. T- 61

Title: Traces of calcium in brine by photometric titration

Summary: Determination of traces of calcium in brine by photometric titration with 1,2-diaminocyclohexanetetraacetic acid using the 610 nm Spectrode.

Sample: Brine

Sample Preparation: none

Instruments and Accessories: 702, 716, 736, 751 or 785 Titrino, 703 Titration Stand, 6.5501.01X Spectrode 610 nm, Metrodata TiNet 2

Analysis: Pipette 50 mL sample into an Erlenmeyer flask, then add 100 μL w(KOH) = 50%, 150 μL triethanolamine and 0.6 mL hydroxynaphthol blue colour indicator (50 mg / 50 mL ultrapure water). Remove the dissolved gases by applying a vacuum for about 1 min, then titrate with 1,2-diaminocyclohexanetetraacetic acid c(DCTA) = 0.062 mmol/L using the MET mode.

Standard addition for linear regression:
Repeat the determination with new samples to which 4 μg , 8 μg or 12 μg Ca^{2+} (standard solution) has been added.

Calculation: a) Direct determination:
1 mL c(DCTA) = 0.062 mmol/L corresponds to 2.48 μg Ca^{2+}

$$\mu\text{g/L Ca}^{2+} = \text{EP1} * \text{C01} * \text{C02} / \text{C00}$$

EP1 = titrant consumption in mL (intersection point TiNet)

C00 = 50 (sample volume in mL)

C01 = 1000 (conversion factor in mL/L)

C02 = 2.48

b) Determination by standard addition (linear regression):
 Calculation by means of linear regression is only possible with
 TiNet 2.

Regression line: $y = a * x + b$

a = slope in mL/ μ g

b = axis intercept in mL

x = addition of Ca^{2+} (standard) in μ g

y = titrant consumption in mL (intersection point TiNet)

$$\mu\text{g/L Ca}^{2+} = b / a * C01 / C00$$

C00 = 50 (sample volume in mL)


C01 = 1000 (conversion factor in mL/L)

Results:

Determination of calcium in 30% KCl brine, prepared from «su-
 rapur» chemicals:

a) EP1 = 0.961 mL \Rightarrow 47.8 $\mu\text{g/L Ca}^{2+}$

b) Linear regression



a = 0.387 mL/ μ g

b = 0.9144 mL

\Rightarrow 47.3 $\mu\text{g/L Ca}^{2+}$