

Ti Application Note No. T- 62

Title:	Analysis of nitrite solutions
Summary:	Determination of nitrite in aqueous solutions by potentiometric back-titration of the added permanganate excess with ammonium iron(II) sulphate using the Pt Titrode.
Sample:	E.g. solution of ammonium nitrite
Sample Preparation:	none
Instruments and Accessories:	702, 716, 736, 751 or 785 Titrino or 726 Titroprocessor, 6.0431.100 Pt Titrode
Analysis:	<p>Blank: Pipette 30.0 mL $c(1/5 \text{ KMnO}_4) = 0.2 \text{ mol/L}$ into a beaker. Add 20 mL $w(\text{H}_2\text{SO}_4) = 25\%$ and titrate with $c[(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2] = 0.2 \text{ mol/L}$.</p> <p>Sample: Pipette 30.0 mL $c(1/5 \text{ KMnO}_4) = 0.2 \text{ mol/L}$ and 20 mL $w(\text{H}_2\text{SO}_4) = 25\%$ into a beaker. Add, e.g., 7.50 mL sample while stirring for 5 min, then titrate with $c[(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2] = 0.2 \text{ mol/L}$.</p>
Calculation:	<p>1 mL $c[(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2] = 0.2 \text{ mol/L}$ corresponds to 6.404 mg NH_4NO_2 4.601 mg NO_2^-</p> <p>$\text{g/L NH}_4\text{NO}_2 = (\text{C01} - \text{EP1}) * \text{C02} / \text{C00}$</p> <p>EP1 = titrant consumption for the sample in mL C00 = 7.50 (sample volume in mL) C01 = titrant consumption for the blank in mL C02 = 6.404</p>
Results:	AVG(7) = $136.95 \pm 0.71 \text{ g/L NH}_4\text{NO}_2$