

Ti Application Note No. T- 64

Title:	Titanium and iron in mixtures
Summary:	Simultaneous determination of titanium and iron by potentiometric titration with potassium dichromate using a platinum electrode. Prior to the determination Ti^{4+} and Fe^{3+} are reduced with Cr^{2+} .
Sample:	Solution containing Ti^{4+} and Fe^{3+}
Sample Preparation:	Pipette a suitable volume of sample solution into a beaker, then add 40 mL acid mixture (containing 75 mL conc. HCl and 168 mL conc. H_2SO_4 per litre). Pass nitrogen through the solution and, while stirring, add about 0.2 g solid $CrCl_2$.
Instruments and Accessories:	702, 716, 736, 751 or 785 Titrino or 726 Titroprocessor, 6.0330.100 Pt electrode, 6.0726.100 Ag/AgCl reference electrode
Analysis:	<p>Pass nitrogen through the solution and titrate with $c(1/6 K_2Cr_2O_7) = 0.1$ mol/L just past the third equivalence point.</p> <p>The first equivalence point of the titration curve corresponds to the Cr^{2+} excess, the second EP corresponds to Ti^{3+} (which is oxidised to Ti^{4+}) and the third EP to Fe^{2+} (which is oxidised to Fe^{3+}).</p>
Calculation:	$g/L Ti = (EP2 - EP1) * C01 * C02 / C00$ $g/L Fe = (EP3 - EP2) * C01 * C03 / C00$ <p>EP1 = titrant consumption to reach the first EP in mL EP2 = titrant consumption to reach the second EP in mL EP3 = titrant consumption to reach the third EP in mL C00 = sample volume in mL C01 = 0.1 («normality» of the titrant in mol/L) C02 = 47.88 [M(Ti) in g/mol] C03 = 55.85 [M(Fe) in g/mol]</p>
Results:	$AVG(3) = 119.93 \pm 0.11$ g/L Ti $AVG(3) = 53.17 \pm 0.17$ g/L Fe