

Ti Application Note No. T- 69

Title:	Iron and nickel in binary mixtures
Summary:	Determination of iron and nickel in binary mixtures by potentiometric titration with EDTA at different pH values using the Cu ISE.
Sample:	Solutions containing Fe^{3+} and Ni^{2+} or dissolved alloys
Sample Preparation:	none
Instruments and Accessories:	702, 716, 736, 751 or 785 Titrino or 726 Titroprocessor, 6.0502.140 Cu ISE, 6.0726.100 Ag/AgCl reference electrode
Analysis:	<p>a) Determination of nickel: Pipette, e.g., 3.0 mL sample solution into a beaker. Add 20 mL dist. water, 10 mL $\text{NH}_3/\text{NH}_4\text{Cl}$ buffer pH = 10 and 1 mL CuEDTA complex solution (0.05 mol/L). Stir for 25 s, then titrate with $c(\text{EDTA}) = 0.1 \text{ mol/L}$.</p> <p>b) Determination of iron: Pipette, e.g., 3.0 mL sample solution into a beaker. Add 20 mL dist. water, 10 mL acetate buffer pH = 4.7 and 10.0 mL $c(\text{EDTA}) = 0.1 \text{ mol/L}$. Stir for 25 s, then titrate back the EDTA excess with $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$.</p>
Calculation:	<p>a) 1 mL $c(\text{EDTA}) = 0.1 \text{ mol/L}$ corresponds to 5.869 mg Ni^{2+}</p> $\text{g/L Ni}^{2+} = \text{EP1} * \text{C01} / \text{C00}$ <p>EP1 = titrant consumption in mL C00 = 3.0 (sample volume in mL) C01 = 5.869</p> <p>b) 1 mL $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ corresponds to 5.585 mg Fe^{3+}</p> $\text{g/L Fe}^{3+} = (\text{C01} - \text{EP1}) * \text{C02} / \text{C00}$ <p>EP1 = titrant consumption in mL C00 = 3.0 (sample volume in mL) C01 = 10.0 [added volume of $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ in mL] C02 = 5.585</p>

Results:	AVG(4) = 9.02 ± 0.05 g/L Fe^{3+} AVG(4) = 3.41 ± 0.06 g/L Ni^{2+}
-----------------	--

Remarks:	In the same way binary mixtures of $\text{Fe}^{3+}/\text{Co}^{2+}$ can be analysed.
-----------------	---