


Carboxyl end groups in polymers

Photometric determination based on ASTM D7409



The carboxyl end groups (CEG) in polymers, such as polyethylene terephthalate (PET), are a measure of the number of unreacted carboxylic acid groups at each end of a polymer chain. The number of CEGs may influence the hydrolysis resistance of geosynthetics, such as geogrids and geotextiles. The lower the CEG value the higher is the hydrolysis resistance of geosynthetics, which in turn increases their stability.

This Application Note describes the photometric titration of carboxyl end groups in PET pellets using the Metrohm Optrode. The acidic end groups of the polymer are titrated with an ethanolic KOH solution using bromophenol blue as indicator.

Method description

Sample

Polyethylene terephthalate (PET) pellets

Sample preparation

Approximately 1.0–1.6 g PET pellets are weighed into a 50 mL Erlenmeyer flask and 20 mL o-cresol is added. It might be that the sample quantity and the procedure have to be adapted to prevent precipitation. The solution is boiled while stirring under reflux until the sample is completely dissolved (approx. 20 minutes).

The solution is cooled down to room temperature and 35 mL chloroform and 0.04 mL bromophenol blue are added.

Configuration

905 Titrand	2.905.0010
800 Dosino	2.800.0010
Dosing unit 5 mL	6.3032.150
801 Stirrer	2.801.0040
Optrode ($\lambda = 610$ nm)	6.1115.000
Reflux condenser	
Erlenmeyer flask 50 mL	

Solutions

o-Cresol p.a.	If crystalline, liquefy carefully in a dry box at 35 °C. CAS 95-48-7
Chloroform	CAS 67-66-3
Titrand	c(KOH) = 0.05 mol/L in ethanol CAS 1310-58-3
Bromophenol blue indicator	Dissolve approx. 0.5 g bromophenol blue in 100 mL methanol CAS 115-39-9
o-Cresol/chloroform mixture	In a ratio of 2:3 (v/v)

Analysis

Blank

The sample preparation is performed as for the PET sample but this time without sample. After the solution has cooled down to room temperature, 0.04 mL bromophenol blue indicator is added and the solution is titrated from yellow to blue with ethanolic KOH.

Sample

To the dissolved sample, 0.04 mL bromophenol blue indicator is added and the solution is titrated from yellow to blue with ethanolic KOH.

Parameters

Titration mode	MET U
Measurement drift	50 mV/min
Min. waiting time	0 s
Max. waiting time	26 s
Volume increment	0.01 mL (for blank) 0.1 (for titer and sample)
Dosing rate	5 mL/min
Pause	40 s
Filling rate	2.5 mL/min
EP criterion	30 mV
EP recognition	greatest
Stirring rate	8

Results

Carboxyl end groups / (mmol/kg) (n = 5)	s(rel) / %
32.11	2.1