


Titration Application Note T-098

Conductometric determination of TBN in petroleum products according to IP 400



This Application Note describes the conductometric determination of the total base number according to IP 400. The investigated sample is hydraulic oil.

Method description

Sample

Hydraulic oil

Sample preparation

No sample preparation required

Configuration

| | |
|--|------------|
| 856 Conductivity module | 2.856.0220 |
| 800 Dosino 2 × | 2.800.0010 |
| Dosing unit 10 mL | 6.3032.210 |
| Dosing unit 50 mL | 6.3032.250 |
| 804 Ti-Stand | 2.804.0010 |
| 802 Rod Stirrer | 2.802.0020 |
| Conductivity cell, $c = 0.1 \text{ cm}^{-1}$ | 6.0916.040 |

Solutions

| | |
|---------|--|
| Titrant | $c(\text{HCl}) = 0.1 \text{ mol/L}$ in 2-propanol. If possible this solution should be bought from a supplier |
| Solvent | $\Phi(\text{toluene}) = 50\%$, $\Phi(2\text{-propanol}) = 49.5\%$ and $\Phi(\text{CO}_2\text{-free H}_2\text{O}) = 0.5\%$ |

Analysis

An appropriate amount of well-mixed sample is weighed into the titration vessel and 75 mL solvent are added. Proceed with the titration using standardized $c(\text{HCl}) = 0.1 \text{ mol/L}$.

After the titration, the burette tip and conductivity measuring cell are rinsed with solvent until no more oil is visible. The conductivity measuring cell is then rinsed with dist. H_2O followed by solvent.

Parameters

| | |
|-------------------|----------|
| Mode | MET Cond |
| Stirrer speed | 8 |
| Pause | 30 s |
| Signal drift | off |
| Min. waiting time | 10 s |
| Max. waiting time | 10 s |
| Volume increment | 0.1 mL |

Calculations

$$\text{TBN} = \frac{(V_{\text{last EP}} - V_{\text{blank}}) \times c_{\text{HCl}} \times f \times M_{\text{A}}}{m_{\text{s}}}$$

| | |
|----------------------|---|
| TBN | Total acid number in mg KOH/g sample |
| $V_{\text{last EP}}$ | Titration consumption in mL to reach the last equivalence point. |
| V_{blank} | Blank value consumption for the used quantity of solvent |
| c_{HCl} | Concentration of titrant in mol/L; here $c(\text{HCl}) = 0.1 \text{ mol/L}$ |
| f | Correction factor (titer), dimensionless |
| M_{A} | Molar mass of KOH; 56.106 g/mol |
| m_{s} | Sample size in g |

Results

| | Sample size/(g) | TBN/ (mg KOH/g sample) |
|--------|-----------------|---------------------------|
| 1 | 3.0378 | 8.206 |
| 2 | 3.0336 | 8.256 |
| 3 | 3.0181 | 8.219 |
| 4 | 3.0088 | 8.193 |
| 5 | 3.0124 | 8.198 |
| 6 | 3.0024 | 8.197 |
| Mean | | 8.212 |
| s(rel) | | 0.29% |

www.metrohm.com

 **Metrohm**