


Titration Application Note T-114

Iodometric determination of ascorbic acid in orange juice



This application note describes the iodometric, bivalent determination of ascorbic acid in orange juice using the Double Pt-sheet electrode

Method description

Sample

Orange juice

Blood orange juice

Sample preparation

No sample preparation is required.

Configuration

907 Titrand	2.907.0010
801 Magnetic Stirrer	2.801.0040
800 Dosino (3)	2.800.0010
5 mL Dosing unit	6.3032.150
10 mL Dosing unit	6.3032.210
20 mL Dosing unit	6.3032.220
Double Pt-sheet electrode	6.0309.100

Solutions

Titrand	$c(I_2) = 0.01 \text{ mol/L}$
Sodium hydroxide	$c(\text{NaOH}) = 1 \text{ mol/L}$
Sulfuric acid	$w(\text{H}_2\text{SO}_4) = 25\%$
Glyoxal solution $w(\text{glyoxal}) = 40\%$	200 mL $w(\text{glyoxal}) = 40\%$ is adjusted with $c(\text{NaOH}) = 1 \text{ mol/L}$ to $\text{pH} = 7.0$. The solution is stored away from light in a dark bottle in a refrigerator.

Analysis

50 mL juice (0.05–0.5 mg ascorbic acid) is pipetted into a titration beaker; 2 mL glyoxal solution is added, and the mixture is stirred briefly and allowed to stand for 5 min. After the addition of 5 mL sulfuric acid $w(\text{H}_2\text{SO}_4) = 25\%$, the sample is titrated with $c(I_2) = 0.01 \text{ mol/L}$ until after the equivalence point.

Parameters

Mode	MET Ipol
Vol. increment	0.1 mL
Signal drift	50 mV/min
Max. waiting time	26 s
I(pol)	1 μA
EP criterion	30 mV

EP recognition	greatest
----------------	----------

Results

Mean result for orange juice (n = 3)

Ascorbic acid / (mg/L)	s(rel)
318.97	0.96%

Mean result for blood orange juice (n = 3)

Ascorbic acid / (mg/L)	s(rel)
470.78	0.79%