


# Bivoltametric determination of ascorbic acid in orange juice using 2,6-dichlorophenol-indophenol



This Application Note describes the bivoltametric determination of ascorbic acid in orange juice using 2,6-dichlorophenolindophenol (DPIP) as titrant and the Double Pt-sheet electrode for indication.

# Method description

## Sample

Orange juice

## Sample preparation

No sample preparation is required.

## Configuration

907 Titrand	2.907.0010
801 Magnetic Stirrer	2.801.0040
800 Dosino (2)	2.800.0010
50 mL Dosing unit	6.3032.250
20 mL Dosing unit	6.3032.220
Double Pt sheet electrode	6.0309.100

## Solutions

Titrand c(DPIP) = 0.25 g/L	50 mg DPIP sodium salt dihydrate is dissolved in about 150 mL hot dist. water (50 – 60 °C) containing 42 mg sodium hydrogen carbonate. The solution is transferred into a 200 mL volumetric flask and made up to the mark with dist. H <sub>2</sub> O. The solution can be stored up to one week in a brown glass bottle in a refrigerator. The titer of this solution has to be determined daily.
Oxalic acid solution w(H <sub>2</sub> C <sub>2</sub> O <sub>4</sub> ) = 2%	28 g oxalic acid dihydrate is weighed into a 1 L volumetric flask and dissolved in dist. H <sub>2</sub> O. The flask is then filled up to the mark with dist. H <sub>2</sub> O.

## Analysis

1 mL sample (0.05–0.5 mg ascorbic acid) is pipetted into a titration beaker. 50 mL oxalic acid solution is added. The solution is then titrated with β(DPIP) = 0.25 g/L until after the equivalence point.

## Parameters

Mode	MET Ipol
Stirring rate	4
Start volume	1 mL
Pause	30 s
Vol. increment	0.05 mL
Signal drift	30 mV/min
Max. waiting time	32 s
I(pol)	1 μA
EP criterion	30 mV
EP recognition	greatest

## Results

Mean result for orange juice (n = 3)

Ascorbic acid / (mg/L)	s(rel)
372.00	0.80%