


Titration Application Note T-138

# Determination of inorganic sulfate in a secondary alkyl sulfonate according to DIN EN 14880



Inorganic sulfate in a secondary alkyl sulfonate raw material is determined according to DIN EN 14880 using the Pb-ISE for indication.

# Method description

## Sample

Secondary alkyl sulfonate raw material, purity approx. 60%

## Sample preparation

Approx. 1.4 g sample is weighed accurately into a 100 mL volumetric flask. The sample is then dissolved in deion. H<sub>2</sub>O. To prevent foam formation, 15 mL isopropanol is added to the solution and the flask is filled up to the mark with deion. H<sub>2</sub>O.

## Configuration

907 Titrand	2.907.0020
801 Magnetic stirrer	2.801.0040
800 Dosino, 3x	2.800.0010
50 mL Dosing unit	6.3032.250
20 mL Dosing unit	6.3032.220
5 mL Dosing unit	6.3032.150
Ion-selective electrode Pb	6.0502.170
LL ISE Reference	6.0750.100

## Solutions

Titrand	$c(\text{Pb}(\text{NO}_3)_2) = 0.01 \text{ mol/L}$ Approx. 3.31 g $\text{Pb}(\text{NO}_3)_2$ is weighed into a 1 L volumetric flask and dissolved in deion. H <sub>2</sub> O. The flask is then filled up to the mark with deion. H <sub>2</sub> O.
Solvent	Isopropanol
Nitric acid	$c(\text{HNO}_3) = 1 \text{ mol/L}$

## Analysis

10 mL sample solution is pipetted into a titration beaker. After the addition of 50 mL isopropanol and 0.5 mL  $c(\text{HNO}_3) = 1 \text{ mol/L}$ , the solution is titrated with  $c(\text{Pb}(\text{NO}_3)_2) = 0.01 \text{ mol/L}$  until after the first equivalence point.

## Parameters

Mode	DET U
Stirring rate	8
Start volume	1.5 mL
Pause	30 s
Signal drift	20 mV/min
Max. waiting time	38 s
Meas. point density	2
Min. increment	10 $\mu\text{L}$
Dosing rate	2 mL/min
EP criterion	5
EP recognition	greatest

## Results

Content of sulfate as sodium sulfate in %

$w(\text{Na}_2\text{SO}_4) / \%$	$s(\text{rel}), n = 5$
4.32	0.35%